TABLE OF CONTENTS

English...	1	German...	5
French...	2	Czech...	6
Spanish...	3	Product Codes...	7
Italian...	4	Glossary of Symbols...	7

INTENDED USE

QC1™ Gram Quality Control Slides are designed for the quality control of the Gram stain technique. They utilize Gram-positive bacteria and Gram-negative bacteria as Gram stain reaction controls. QC1 Gram Slides are intended to be used to validate each batch of slides stained by the conventional Gram stain technique.

SUMMARY

Various clinical laboratory regulatory organizations and Good Laboratory Practices require quality control programs. The performance of equipment, reagents, methods, and techniques must be monitored and documented. The QC1 Gram Slide is a control slide for use with the Gram Stain test, eliminating the need to maintain in-house stock cultures for preparing slides.

FOR IN VITRO DIAGNOSTIC USE ONLY

REAGENT COMPONENTS

QC1 Gram Slides contain substrates (bacterial cell suspensions) at optimal concentrations (CFU/ml) necessary for a positive and a negative Gram stain reaction. The positive control, *Staphylococcus aureus*, displays a Gram-positive staining reaction and the negative control, *Escherichia coli*, displays a Gram-negative staining reaction.

PRECAUTIONS

Handle all slides by the edges. Do not touch the surface of the slide as this may damage the substrates. The control slides contain bacterial concentrations that are air-dried and heated to adhere to the slide surface. Routine bacteriological safety controls should be used when handling the bacterial slides.

STABILITY AND STORAGE

The slides are stable to the stated expiration date when stored at room temperature (15–30°C).

PROCEDURE

Materials Provided: QC1 Gram Slides.

Materials Not Provided: Gram staining reagents, slides, microscope with oil immersion lens, and absorbent material.

Method of Staining: The substrates on the QC1 Gram Slides are pre-fixed and do not need to be fixed by heat or alcohols prior to staining. If using QC1 Gram Slides that have an area for a patient specimen, it is acceptable to fix [heat] the patient area[s] of the slide by using appropriate temperatures. Add the QC1 Gram Slide at the first stain step of the staining procedure. Stain following the procedure described by the manufacture of the stain being used or consult appropriate references.

EXPECTED RESULTS

Retention of the purple mordant primary stain indicates a Gram-positive organism. Bacterial cells that decolorize and stain pink/red with counterstain are Gram-negative organisms.

LIMITATIONS OF PROCEDURES

For patient slides, a proper heat fixing technique must be used to prevent change of the cell wall and the chemical structure of the cell wall. Care must be taken not to heat the slide for too long or too hot. A microorganism that is physically disrupted by excess heating will not react to Gram staining as expected.

BIBLIOGRAPHY

1. Todd & Sanford, Clinical Diagnosis By Laboratory Methods, 14th Edition.
2. Stain Technology, Vol. VIII, No. 2, January 1932.
3. Manual For Microbial Methods, NY McGraw-Hill, 1957, pp. 15–18.
4. Diagnostic Microbiology, 2nd Edition, 1966, p. 320.

CONTACT

For Technical Assistance email Technical@AlphaTecSystems.com and for Customer Service, email Sales@AlphaTecSystems.com or call [+1] 800.221.6058 or [+1] 360.260.2779 between 8 am and 4 pm Monday through Friday, Pacific Time.

WARRANTY

CalibreScientific US, Inc. Inc. warrants this product to perform as described in the labeling and literature supplied. CalibreScientific US, Inc. disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall CalibreScientific US, Inc. be liable for any consequential damages arising out of aforesaid express warranty.

TRADEMARKS

QC1™ is a trademark of CalibreScientific AMER IP LLC., 6201 Trust Dr, Holland, OH 43528.

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QC1™ je ochranná známka společnosti CalibreScientific AMER IP LLC., 6201 Trust Dr, Holland, OH 43528.
Manufactured by CalibreScientific US, Inc.
MDSS GmbH