## TABLE OF CONTENTS

English...	1	Portuguese...	6
French...	2	Russian...	7
Spanish...	3	Chinese...	8
Italian...	4	Product Codes...	9
German...	5	Glossary of Symbols...	9

## INTENDED USE

Alpha-Tec QC1™ AFB Quality Control Slides are designed for the quality control of the acid-fast bacilli stain techniques. They utilize an acid-fast positive bacteria and an acid-fast negative bacteria as acid-fast stain reaction controls. QC1 AFB Slides are intended to be used daily to validate each batch of slides stained and are compatible with the various acid-fast stain techniques (Ziehl-Neelsen, Kinyoun Carbolfuchsin, and fluorescent procedures).

## SUMMARY

Quality control programs are required by various regulatory organizations and Good Laboratory Practices to monitor and document the performance of equipment, reagents, methods, and techniques. QC1 AFB Slides are for use with the acid-fast stain test, and eliminate the need to maintain in-house stock cultures for preparing slides.

## FOR IN VITRO DIAGNOSTIC USE ONLY

## TEST COMPONENTS

QC1 AFB Slides utilize substrates that consist of bacterial cell suspensions at optimal concentrations for microscopically viewing the correct staining results. The positive well (+) contains *Mycobacterium scrofulaceum* to produce the positive acid-fast stain reaction. The negative well (-) contains *Escherichia coli*, which results in the negative acid-fast stain reaction.

## STABILITY AND STORAGE

The substrates on QC1 AFB Slides are stable to the stated expiration date when stored at room temperature (15–30°C).

## PRECAUTIONS

Handle all slides by the edges. Do not touch the surface of the slide as this may damage the substrates. Control slides contain bacterial concentrations that are air-dried and heated to adhere to the slide surface. Routine bacteriological safety controls should be used when handling the bacterial slides.

## PROCEDURE

**Materials Provided:** QC1 AFB Slides

**Materials Not Provided:** Staining reagents and microscope with oil immersion lens.

**Method of Staining:** Alpha-Tec QC1 AFB Slides are pre-fixed and should **not** be fixed by heat or alcohols prior to staining. Add QC1 AFB Slides to the first stain step of the staining procedure. Stain QC1 AFB Slides according to the procedure described by the manufacturer of the stain, or consult appropriate references. **Note:** Alpha-Tec offers a complete line of acid-fast stains. Contact Customer Service for additional information.

**Technical Note:** In order to prevent any potential cross-contamination of the control wells between staining steps, rinse the slide by tilting it onto its side not on its end.

## RESULTS

Refer to the instruction sheet provided with the stain of choice or consult appropriate references.

## LIMITATIONS OF THE PROCEDURE

A positive staining reaction with an acid-fast stain is presumptive information indicating the presence of mycobacteria in the specimen. For positive identification, culture methods must be employed.

## BIBLIOGRAPHY

1. Procedures for the isolation and identification of mycobacteria. U.S. Dept. of Health, Education and Welfare (CDC). 1975. #79–8230.
2. Laboratory Manual for Acid-Fast Microscopy. 2nd Ed. U.S. Dept. of Health, Education and Welfare (CDC). 1976.
3. Lennette, E.H., Manual of Clinical Microbiology. 3rd Ed. 1980, pp. 150–179.

## CONTACT

Alpha-Tec Systems, Inc. offers a complete line of digestion and decontamination systems, stains and reagents for pre-analytical specimen processing. For Technical Assistance email [Technical@AlphaTecSystems.com](mailto:Technical@AlphaTecSystems.com) and for Customer Service, email [Sales@AlphaTecSystems.com](mailto:Sales@AlphaTecSystems.com) or call [+1] 800.221.6058 or [+1] 360.260.2779 between 8 am and 4 pm Monday through Friday, Pacific Time.

## WARRANTY

Alpha-Tec Systems, Inc. warrants this product to perform as described in the labeling and literature supplied. Alpha-Tec Systems, Inc. disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall Alpha-Tec Systems, Inc. be liable for any consequential damages arising out of aforesaid express warranty.

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QC1™ is a trademark of Alpha-Tec Systems, Inc., 1311 SE Cardinal Court, Suite 170, Vancouver, WA 98683 USA.
Manufactured by Alpha-Tec Systems, Inc.
MDSS GmbH