TABLE OF CONTENTS

English...	1	Product Codes...	4
Spanish...	3	Glossary of Symbols...	4

INTENDED USE

QC1™ PROTO-FIX® Slides were designed for the quality control of a Modified Wheatley's Trichrome Stain technique utilizing PROTO-FIX fixative. The slides utilize a controlled parasite (*Giardia lamblia*) fixed in PROTO-FIX yielding a positive staining reaction. QC1 PROTO-FIX Slides are intended to validate the quality of the staining reagents and the technique of the procedure.

SUMMARY

The identification of the intestinal protozoa is based upon unstained and stained morphology, as observed in several types of staining preparations. The QC1 PROTO-FIX Slides, used in conjunction with the Modified Wheatley's Trichrome Stain procedure (#0003351), provide convenient, rapid and effective means for the quality control of this stain procedure.

FOR IN VITRO DIAGNOSTIC USE ONLY

PRECAUTIONS

Handle the substrate slides by the edges. Do not touch the surface of the slide as this may damage the substrate. Observe all safety precautions for handling stool specimens.

STABILITY AND STORAGE

The substrate slides are stable to the stated expiration date when stored at the required temperature. Store at room temperature (15-30°C). Avoid excessive heat and sunlight.

USER QUALITY CONTROL

Using good laboratory practices, new lot numbers of quality control slides should be stained and compared against a previous lot number of the same product number. Also, all changes in staining reagents should be validated against a known quality control slide.

SPECIMEN COLLECTION & PREPARATION

Handle the slides from the edges only. The substrate, which contains the quality control organisms, appears clear and colorless. Touching, rubbing, or smearing of the top surface of the slide could result in removal of the quality control substrate resulting in an unacceptable stained result.

PROCEDURE

Materials Provided: QC1 PROTO-FIX Slides

Materials Not Provided: Coverslips, Microscope, staining reagents of choice.

Modified Wheatley's Trichrome Stain Procedure

Reagent	Timing
70% Ethanol (#0003359)	1.5 minutes. Remove and drain off excess liquid by touching edge of slide to absorbent material.
70% Ethanol	1.5 minutes. Remove and drain off excess liquid by touching edge of slide to absorbent material.
Wheatley's Trichrome Stain* (#0003351)	13 minutes. Remove and drain off excess liquid by touching edge of slide to absorbent material.
90% Acid-Ethanol (#0003350)	5 seconds. Remove and drain off excess liquid by touching edge of slide to absorbent material.
100% Ethanol** (0003303)	10 seconds. Remove and drain off excess liquid by touching edge of slide to absorbent material.

Reagent	Timing
100% Ethanol**	1 minute. Remove and drain off excess liquid by touching edge of slide to absorbent material.
100% Ethanol**	1 minute. Remove and drain off excess liquid by touching edge of slide to absorbent material.
Xylene (#0003342) or PRO-CLEAR™ (#0003336)	3 minutes. Remove and drain off excess liquid by touching edge of slide to absorbent material and coverslip with mounting fluid. Read under oil immersion.

NOTE: For best results, scan the QC1 PROTO-FIX Slides first under low power, under high power, and then oil immersion at 100X.

***NOTE:** Staining time varies depending on the intensity desired for the final stain.

****NOTE:** 100% Ethanol or ethanol based reagent alcohol (or equivalent).

EXPECTED RESULTS

These slides are made with cultured *Giardia lamblia* trophozoites and will exhibit slightly different colors than *Giardia lamblia* fixed in PROTO-FIX sourced from clinical fecal samples. The cytoplasm of the trophozoites will be bluish-purple (blue-violet). The chromatin will stain purplish-red (red-violet). Background material and artifacts are not usually seen.

LIMITATIONS OF PROCEDURES

- QC1 PROTO-FIX Slides are prepared with cultured *Giardia lamblia* trophozoites. Therefore, there will be no fecal debris or other artifacts in the field that are usually seen in clinical specimens. The stain results on the control slide will have less intensity than stained *Giardia lamblia* from clinical intestinal specimens.
- The xylene substitute AmeriClear® does not yield appropriate stain results and should not be used with PROTO-FIX.

BIBLIOGRAPHY

- Allen, K, Frankel, JW, and Kwa, B. Comparison of CONSED® RED and formalin-ethyl acetate methods for concentrating intestinal parasites and eggs. Southeastern Branch ASM. 1997 meeting.
- Amin, O. Evaluation of a new system for the fixation, concentration, and staining of intestinal parasites in fecal specimens, with critical observations on the Trichrome stain. Journal of Microbiological Methods. 2000. 38:127-132.
- Erdman, D. Clinical comparison of ethyl acetate and diethyl ether in the formalin-ether sedimentation technique. J Clin Micro. 1981. 14:483-485.
- Faust, EC et al. A critical study of clinical laboratory techniques for the diagnosis of protozoan cysts and helminth eggs in feces. AM J Trop. Hyg Med. 1938. 18:169-183.
- Garcia, LS and Ash, L. Diagnostic Parasitology. C.V. Mosby Co. St. Louis, MO. 1979. pp.9-24.
- Jensen, B. et al. Comparison of polyvinyl alcohol fixative with three less hazardous fixatives for detection and identification of intestinal parasites. J Clin Micro. 2000. 38(4): 1592-1598.
- Markell, E and Voge, M. Medical Parasitology. W.B. Saunders Co. Philadelphia, PA. 1981. pp.317-355.
- Melvin, DM and Brooke, MM. Laboratory procedures for the diagnosis of intestinal parasites. 3rd Edition. US Department of Health and Human Services Pub. No. (CDC). 1982. 82-8282.
- Price, DL. Procedure manual for the diagnosis of intestinal parasites. CRC Press, Inc. Cleveland, OH. 1994.
- Sholten, TH and Yang, J. Evaluation of unpreserved and preserved stools for the detection and identification of intestinal parasites. AM J Clin Pathol. 1974. 62:563-567.
- Williams, RW and Dennis, M. Technical Performance of a new fecal preservative, (PROTO-FIX) compared with 10% formalin for

ELISA Assays for *Giardia lamblia* and *Cryptosporidium parvum* in fecal specimens. ASM Meeting, Atlanta, GA. 1998.

12. Young, K, Bullock, S, Melvin, D, and Spruill, C. Ethyl acetate as a substitute for diethyl ether in the formalin-ether sedimentation technique. J Clin Micro. 1979. 10:852-853.

CONTACT

Alpha-Tec Systems, Inc. offers a complete line of reagents, stains, and QC1™ Quality Control Slides for AFB, Parasitology, Bacteriology, and Mycology processing, as well as O&P collection systems and concentration devices for Parasitology. For Technical Assistance, email Technical@AlphaTecSystems.com, and for Customer Service, email Sales@AlphaTecSystems.com, or call either [+1] 800.221.6058 (USA) or [+1] 360.260.2779 between 8AM and 4PM Monday through Friday, Pacific Time.

WARRANTY

This product is warranted by Alpha-Tec Systems, Inc. to perform as described in the labeling and literature supplied. Alpha-Tec Systems, Inc. disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall Alpha-Tec Systems, Inc. be liable for any consequential damages arising out of aforesaid express warranty.

TRADEMARKS

CONSED®, PRO-CLEAR™, PROTO-FIX® and QC1™ are trademarks of Alpha-Tec Systems, Inc., 1311 SE Cardinal Court, Suite 170, Vancouver, WA 98683 USA.

AmeriClear® is a trademark of Cardinal Health, McGaw Park, IL, USA.

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- Allen, K., Frankel, J.W. and Kwa, B., 1997. "Comparison of CONSED® Red and Formalin-ethyl Acetate Methods for Concentrating Intestinal Parasites and Eggs". Southeastern Branch, ASM, 1997 meeting.
- Amin, O., 2000. "Evaluation of a new system for the fixation, concentration, and staining of intestinal parasites in fecal specimens, with critical observations on the Trichrome stain." Journal of Microbiological Methods 38:127-132.

3. Erdman, Dean, 1981. "Clinical Comparison of Ethyl Acetate and Diethyl Ether in the Formalin-ether Sedimentation Technique." *J Clin Micro* 14:483-485.
4. Faust, E.C. et al, 1938. "A Critical Study of Clinical Laboratory Techniques for the Diagnosis of Protozoan Cysts and Helminth Eggs in Feces." *AM J Trop. Hyg. Med* 18:169-183.
5. Garcia, L.S. and Lawrence Ash, 1979. *Diagnostic Parasitology*, C.V. Mosby Co., St. Louis, MO, pp.9-24.
6. Jensen, B., et al, 2000. "Comparison of Polyvinyl Alcohol Fixative with Three Less Hazardous Fixatives for Detection and Identification of Intestinal Parasites." *J Clin Micro* 38(4): 1592-1598.
7. Markell, Edward and Marietta Voge, 1981. *Medical Parasitology*. W.B. Saunders Co., Philadelphia, PA, pp.317-355.
8. Melvin, D.M. and M.M. Brooke, 1982. "Laboratory Procedures for the Diagnosis of Intestinal Parasites." 3rd Edition. US Department of Health and Human Services Pub. No. (CDC) 82-8282, Atlanta, GA.
9. Price, D.L. 1994. *Procedure Manual for the Diagnosis of Intestinal Parasites*. CRC Press, Inc., Cleveland, OH.
10. Sholten, T.H. and J. Yang, 1974. "Evaluation of Unpreserved and Preserved Stools for the Detection and Identification of Intestinal Parasites." *AM J Clin Pathol* 62:563-567.
11. Williams, R.W. and Dennis, M., 1998. Technical Performance of a New Fecal Preservative, (PROTO-FIX) Compared with 10% Formalin for ELISA Assays for *Giardia lamblia* and *Cryptosporidium parvum* in Fecal Specimens. ASM Meeting, Atlanta, GA.
12. Young, K., S. Bullock, D. Melvin and C. Spruill, 1979. "Ethyl Acetate as a Substitute for Diethyl Ether in the Formalin-ether Sedimentation Technique." *J Clin Micro* 10:852-853.
13. Unpublished data on file.

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PRODUCT CODES

0003255 QC1 PROTO-FIX Slides, 20/Box

Manufactured by Alpha-Tec Systems, Inc.
1311 SE Cardinal Court, Suite 170
Vancouver, WA 98683 USA

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MDSS GmbH