**TABLE OF CONTENTS**

English...	1	Product Codes...	3
French...	2	Glossary of Symbols...	3

**INTENDED USE**

QC1™ Cryptosporidium Quality Control Slides are designed for the quality control of the staining technique and to serve as a comparative aid in the identification of *Cryptosporidium* species. The slides utilize *Cryptosporidium* oocysts as a positive control, yielding a positive identification and staining reaction. A negative staining control is also included.

**SUMMARY**

The identification of *Cryptosporidium* spp. is based upon stained morphology as observed in several types of staining preparations. A number of simple techniques have been developed that are satisfactory for diagnostic work.

**FOR IN VITRO DIAGNOSTIC USE ONLY****PRECAUTIONS**

**CAUTION:** Handle all substrate slides by the edges. Do not touch the surface of the slides as this may damage the substrates. The stain control slides contain 10% formalin-fixed concentrations that are air-dried and heat-fixed. Routine safety controls should be used when handling the slides.

**STABILITY AND STORAGE**

The substrate slides are stable to the stated expiration date when stored at room temperature (15-30°C).

**USER QUALITY CONTROL**

QC1 Cryptosporidium Slides consist of a positive "well" containing optimal concentrations of 10% formalin-fixed *Cryptosporidium* oocysts. This results in a positive staining reaction for the particular stain used and a morphologically correct organism set against a background light in debris for superior differentiation of the organism. A negative control is included to confirm a negative staining reaction.

**PROCEDURE**

**Materials Provided:** QC1 Cryptosporidium Slides with two specimen wells/slide.

**Materials Not Provided:** Staining reagents, microscope with oil immersion lens.

**Method of Staining:**

The Alpha-Tec QC1 Cryptosporidium Slides are pre-fixed and should not be fixed by heat or alcohols prior to staining. Add the QC1 Cryptosporidium Slide at the first stain step of the staining procedure. Fix dried patient smears for 9-10 seconds in absolute methanol. Do not rinse with water. Allow to air-dry.

1. Stain using method of choice.
2. Coverslip using mounting media.
3. Observe using high-dry objective to screen. Switch to oil immersion for confirmation.

**EXPECTED RESULTS****KINYOUN CARBOLFUCHSIN ACID-FAST STAIN RESULTS:**

**Positive Well:** *Cryptosporidium* oocysts yield a brilliant pink to red color against a blue counterstained background, showing typical internal vacuole and material clumped to one side of the 4-6 micron cyst. Empty cyst walls (ghost cysts) and damaged oocysts stain poorly.

**Negative Well:** Substrate will stain blue.

**LIMITATIONS OF PROCEDURES**

Acid-fast staining methods are usually trouble-free procedures when followed as directed. If problems do occur, they will usually be the result of the specimen not being completely dried or diluting of reagents from carryover of the rinsing step to the decolorizer / counterstain step. As a result, screening of the specimen using a high-dry objective will show "blurred" morphology and stained oocysts will give red-pink overcasts only.

**BIBLIOGRAPHY**

1. Bronsdon, MA. Rapid dimethyl sulfoxide-modified acid-fast stain of cryptosporidium oocysts in stool specimens. *Journal of Clinical Microbiology*. 1984. 19:952-953.
2. Current, WL. Human enteric coccidia: cryptosporidium. *Clinical Microbiology Newsletter*. 1985. 7:167-170.
3. Gardner, JA, et al. Comparison of three staining methods for the identification of cryptosporidium in stool specimens." Presented 1985 Annual ASM Convention, (Abstract C20).
4. Price, DL. Procedure manual for the diagnosis of intestinal parasites. CRC Press. ISBN 0-8493-8654-3. 1994. pg. 3.

**CONTACT**

Alpha-Tec Systems, Inc. offers a complete line of reagents, stains, and QC1™ Quality Control Slides for AFB, Parasitology, Bacteriology, and Mycology processing, as well as O&P collection systems and concentration devices for Parasitology. For Technical Assistance, email [Technical@AlphaTecSystems.com](mailto:Technical@AlphaTecSystems.com), and for Customer Service, email [Sales@AlphaTecSystems.com](mailto:Sales@AlphaTecSystems.com), or call either [+1] 800.221.6058 (USA) or [+1] 360.260.2779 between 8AM and 4PM Monday through Friday, Pacific Time.

**WARRANTY**

This product is warranted by Alpha-Tec Systems, Inc. to perform as described in the labeling and literature supplied. Alpha-Tec Systems, Inc. disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall Alpha-Tec Systems, Inc. be liable for any consequential damages arising out of aforesaid express warranty.

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## BIBLIOGRAPHIE

1. Bronsdon, MA. Rapid dimethyl sulfoxide-modified acid-fast stain of cryptosporidium oocysts in stool specimens. Journal of Clinical Microbiology. 1984. 19:952-953.
2. Current, WL. Human enteric coccidia: cryptosporidium. Clinical Microbiology Newsletter. 1985. 7:167-170.
3. Gardner, JA, et al. Comparison of three staining methods for the identification of cryptosporidium in stool specimens." Presented 1985 Annual ASM Convention, (Abstract C20).
4. Price, DL. Procedure manual for the diagnosis of intestinal parasites. CRC Press. ISBN 0-8493-8654-3. 1994. pg. 3.

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**PRODUCT CODES**

0003256 QC1™ Cryptosporidium Slides, 20/Box

Manufactured by Alpha-Tec Systems, Inc.  
1311 SE Cardinal Court, Suite 170  
Vancouver, WA 98683 USA

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QC1™ is a trademark of Alpha-Tec Systems, Inc., 1311 SE Cardinal Court, Suite 170, Vancouver, WA 98683 USA.
QC1™ est une marque d'Alpha-Tec Systems, Inc. 1311 SE Cardinal Court, Suite 170, Vancouver, WA 98683 USA.
MDSS GmbH