Giemsa (Malaria) Stain Methylene Blue Phosphate Pre-Stain Giemsa (Malaria) Stain Buffer

TABLE OF CONTENTS

English...	1	Product Codes...	5
Spanish...	3	Glossary of Symbols...	5

Giemsa (Malaria) Stain & Giemsa stain reagents are intended for the microscopic examination of blood smears to establish the presence of Plasmodium parasites. Persons suspected of having malaria and showing no microscopic presence of Plasmodium parasites should have the blood smears repeated every 12-24 hours for three consecutive days. If the smears remain negative, the diagnosis of malaria is unlikely.

SUMMARY

To establish the presence of malaria (Plasmodium parasites), a blood smear must be prepared and stained with Giemsa stain and examined microscopically. Wright stain will not reliably demonstrate the Plasmodium parasites. Thick smears are more sensitive in detecting the Plasmodium parasites due to the concentration of the blood, allowing a greater volume of blood being examined. Plasmodium parasites are always intracellular and demonstrate (if stained correctly) blue cytoplasm with red chromatin.

Giemsa stain is used to differentiate nuclear and/or cytoplasmic morphology of platelets, RBCs, WBCs, and parasites. The most dependable stain for blood parasites, particularly in thick films, is Giemsa stain containing Azure B. The Giemsa stock stain must be diluted for use with water buffered to a pH of 7.0-7.2. The Giemsa Stain stock solution must be protected from moisture because the staining reaction is oxidative. Therefore, the oxygen in water will initiate the reaction and ruin the stock stain. The aqueous working dilution of stain is good only for 1 day.

FOR IN VITRO DIAGNOSTIC USE ONLY

PRECAUTIONS

For safety considerations, analysis of samples should be done within a certified Class II biological safety cabinet (BSC). Any laboratory procedure involving infectious materials or cultures requires standard personal protective safety equipment, such as latex gloves and laboratory coats (or disposable coats), and eye protection. Decontamination of blood can be accomplished with commercially available bleach solutions, diluted to appropriate working solutions. Countertops, work surfaces, and instruments can be wiped with the bleach solution. Lab supplies (pipettes, loops, microscope slide, etc.) should be immersed in a decontamination solution and then autoclaved.

STABILITY AND STORAGE

The Giemsa (Malaria) Stain, Giemsa (Malaria) Stain Buffer, and Methylene Blue Phosphate Pre-Stain are stable to the listed expiration date on the reagent bottles. Store the Giemsa (Malaria) Stain and Giemsa (Malaria) Stain Buffer at room temperature away from direct light. Store the Methylene Blue Phosphate Pre-Stain at room temperature, however, after opening store at 2-8°C away from direct light.

USER QUALITY CONTROL

- A. The stock Giemsa (Malaria) Stain and Methylene Blue Pre-Stain are blue; the Giemsa (Malaria) Stain Buffer solutions are colorless. All the above referenced solutions/stains should be clear, with no visible contamination. **NOTE:** The Methylene Blue Phosphate Pre-Stain contains no preservatives and should be stored at 2-8°C after opening and checked regularly for clarity.
- B. Check the Giemsa stain reagents, including the pH of the Giemsa (Malaria) Stain Buffer, before each use. If Triton X-100™ has been added to the Giemsa (Malaria) Stain Buffer, do not use a colorimetric method to determine the pH, because the Triton X-100 interferes with the color indicators. Use a pH meter to test any Giemsa (Malaria) Stain Buffer that may contain Triton X-100. The Giemsa (Malaria) Stain Buffer is usable if the pH is within 7.0-7.2.

- C. Prepare and stain films from "normal" blood and microscopically evaluate the staining reactions of the RBCs, platelets, and WBCs.
 1. Macroscopically, blood films appear purplish. If blue, the buffer was too alkaline; if pink to red, the buffer was too acidic.
 2. Microscopically, RBCs appear pinkish-gray, platelets appear dark pink, and WBCs have purple-blue nuclei and lighter cytoplasm. Eosinophilic granules are bright purple-red, and neutrophilic granules are purple. Basophilic stippling within uninfected RBCs is blue.
 3. Slight variations may appear in the colors described above depending on the batch of stain used and the character of the blood itself, but if the various morphological structures are distinct, the stain is satisfactory.

SPECIMEN PREPARATION

Two types of blood films are used. The thick film is most valuable in diagnostic work since it enables the observation of a larger quantity of blood per microscope field. The thin film is less useful in diagnostics but is superior for the observation of morphological details.

PROCEDURE

Materials Provided: Stock Giemsa (Malaria) Stain, Methylene Blue Phosphate Pre-Stain, Giemsa (Malaria) Stain Buffer.

Materials Not Provided:

Coplin jars, pipettes, graduated cylinders, positive QC slides, microscope slides, methyl alcohol, microscope, pH meter, Triton X-100 (optional).

NOTE: Wear gloves when performing any of these procedures.

A. Thin blood films (only)

1. Fix air-dried film in acetone-free absolute methyl alcohol by dipping the film briefly (two dips) in a Coplin jar.
2. Remove, stand slide on end (on a paper towel), and allow to air-dry.
3. Slides should be stained in a stain jar/Coplin jar. Slides stained for malarial parasites **should not** be stained face-up.
4. Stain with diluted (5%) Giemsa (Malaria) Stain for 20 minutes. For a 5% dilution, add 2 ml of stock Giemsa (Malaria) Stain to 38 ml Giemsa (Malaria) Stain Buffer. This diluted stain solution may be added to a Coplin jar. **NOTE:** The addition of 0.01% Triton X-100 is optional. **NOTE:** The Giemsa working solution should be used only on the day it is prepared. **NOTE:** Giemsa stain will not overstain blood smears. The staining time depends somewhat on the strength of the staining solution, but stronger solutions do not decrease time proportionally. Regardless of the strength of the working solution, the minimum staining time is at least 10 minutes. Dilute stains and longer staining times usually give more detailed staining of blood elements and blood parasites.
5. Wash by briefly dipping the slide in and out of a Coplin jar filled with Giemsa (Malaria) Stain Buffer (one or two dips). **NOTE:** Excessive washing will decolorize the film.
6. Let the slide air-dry in a vertical position. Examine when slides are completely dry.

B. Thick blood films (only)

- NOTE:** Thick blood films may be pre-stained in Methylene Blue Phosphate Pre-Stain to obtain better quality results. However, immediate Giemsa staining may be done without the Pre-Stain. If the blood films are not going to be stained immediately or if they are going to be sent to another laboratory for staining, pre-staining is essential.
1. Allow film to air-dry thoroughly for several hours or overnight. Do not dry films in an incubator or by heat, because this will fix the blood and interfere with the lysing of the RBCs. Do not fix in methyl alcohol. **NOTE:** If a rapid diagnosis of malaria is needed, thick films can be made slightly thinner than usual, allowed to dry for 1 hour, and then stained.

Giemsa (Malaria) Stain Methylene Blue Phosphate Pre-Stain Giemsa (Malaria) Stain Buffer

2. Dip the slide into the Methylene Blue Phosphate Pre-stain for three (3) seconds. Remove and touch edge of slide to a paper towel to remove excess stain.
 3. Dip in Giemsa (Malaria) Stain Buffer for three (3) seconds. Remove and stand slide on end to dry. (The thick films will partially dehemoglobinize as they stand. The films do not need to dry before staining but can dry if staining setup is not ready.)
 4. Stain the slide with diluted (2%) Giemsa (Malaria) Stain for 50 minutes. For a 2% dilution, add 2 ml of stock Giemsa (Malaria) Stain to 98 ml of Giemsa (Malaria) Stain Buffer, or adjust the resulting diluted volume as appropriate for your laboratory's use. **NOTE:** The addition of 0.01% Triton X-100 is optional.
 5. Rinse the blood films carefully so as to not wash off the thick film, first by dipping them in the Giemsa (Malaria) Stain Buffer (pH 7.0), and then by placing them into another Coplin jar containing the Stain Buffer (pH 7.0) for 3-4 minutes to complete the rinsing step.
 6. Let air-dry in a vertical position. Examine when slides are completely dry.
- C. **Combination thin and thick blood films**
1. Allow the thick film to air-dry thoroughly.
 2. Fix the thin film by placing only the thin film in methyl alcohol (two dips). Be sure not to get the alcohol or its fumes on the thick film by slightly tilting the slide.
 3. Let air-dry in a vertical position with the *thick film up*. Be sure slide is thoroughly dry before staining. Introducing even a minute amount of methyl alcohol into the stain dilution will interfere with the lysing of the RBCs in the thick films.
 4. Stain the entire slide with diluted (2%) Giemsa (Malaria) Stain for 50 minutes. Place the slide in the stain, *thick film down* to prevent the debris caused by dehemoglobinization from falling onto the thin film. For a 2% dilution, add 2 ml of stock Giemsa (Malaria) Stain to 98 ml of Giemsa (Malaria) Stain Buffer. (**NOTE:** The addition of 0.01% Triton X-100 is optional).
 5. Rinse the thin film by briefly dipping the film in and out of a Coplin jar of Giemsa (Malaria) Stain Buffer (one or two dips). Wash the thick film for 3 to 5 minutes. Be sure that the thick film is immersed but *do not allow the water to cover any part of the thin film*.
 6. Let air-dry in a vertical position with the *thick film down*.

PROCEDURE NOTES

- A. Blood films prepared from venipuncture blood when an anticoagulant is used must be prepared within 1 hour of collection. Otherwise, certain morphological characteristics of both parasites and infected RBCs may be atypical. Also, thick blood films may wash off the slide during staining procedure.
- B. The correct pH for all solutions is important. Solutions with the incorrect pH will prevent certain morphological characteristics (stippling) from being visible and will not give typical nuclear and cytoplasmic colors on the stained film.
- C. Utilize a stain a QC slide per your laboratory's protocol.

RESULTS

- A. If *Plasmodium* organisms are present, the cytoplasm will stain blue and the nuclear material stains red to purple-red.
- B. Schüffner's stippling and other inclusions in the RBCs infected by *Plasmodium spp.* stain red.
- C. Nuclear and cytoplasmic colors that are seen in the malarial parasites will also be seen in the trypanosomes and any intracellular leishmaniae that are present.
- D. The sheath of microfilariae may or may not stain with Giemsa, while the body will usually appear blue to purple.

REPORTING RESULTS

- A. Report any parasite, including the stage(s) seen. Do not use abbreviations.
- B. Any laboratory providing malaria diagnoses should be able to identify *Plasmodium vivax* and *Plasmodium ovale*, even in the absence of Schüffner's stippling.

LIMITATIONS OF THE PROCEDURE

- A. Finding no parasites in one set of blood films does not rule out a parasitic infection.
- B. Examine a minimum of 300 oil immersion (1,000x) fields before reporting no parasites found.
- C. Examine the entire smear under low power (100x) for the presence of microfilariae. Remember that the sheath may not be visible.

BIBLIOGRAPHY

1. Balows, A., et al., 1991. Manual of Clinical Microbiology. ASM, Washington, D.C., Fifth Ed.
2. Brooke, M.M., and A.W. Donaldson, 1950. Use of a surface-active agent to prevent transfer of malarial parasites between blood films during mass staining procedures. *J. Parasitol.* 36:84.
3. Melvin, D. M., and M.M. Brooke, 1955. Triton X-100™ in Giemsa staining of blood parasites. *Stain Technol.* 30:269-275.
4. U.S. Dept. of Health and Human Services, Centers for Disease Control and Prevention.
5. Wilcox, A., 1960. *Manual for the Microscopical Diagnosis of Malaria in Man*. U.S. Public Health Service publication no. 796. U.S. Government Printing Office, Washington, D.C. (Out of print.)

CONTACT

CalibreScientific US, Inc. offers a complete line of reagents, stains, and QC1™ Quality Control Slides for AFB, Parasitology, Bacteriology, and Mycology processing, as well as O&P collection systems and concentration devices for Parasitology. For Technical Assistance, email Technical@AlphaTecSystems.com, and for Customer Service, email Sales@AlphaTecSystems.com, or call either [+1] 800.221.6058 (USA) or [+1] 360.260.2779 between 8AM and 4PM Monday through Friday, Pacific Time.

WARRANTY

This product is warranted by CalibreScientific US, Inc. to perform as described in the labeling and literature supplied. CalibreScientific US, Inc. disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall CalibreScientific US, Inc. be liable for any consequential damages arising out of aforesaid express warranty.

TRADEMARKS

QC1™ is a trademark of CalibreScientific AMER IP LLC., 6201 Trust Dr, Holland, OH 43528..

Triton X-100™ is a trademark of Union Carbide Chemicals & Plastics Company, Inc., 39 Old Ridgebury Road, Danbury, CT 06817.

Giemsa (Malaria) Stain Methylene Blue Phosphate Pre-Stain Giemsa (Malaria) Stain Buffer

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Giemsa (Malaria) Stain Methylene Blue Phosphate Pre-Stain Giemsa (Malaria) Stain Buffer

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1. Balows, A., et al., 1991. Manual of Clinical Microbiology. ASM, Washington, D.C., Quinta Edición.
2. Brooke, M.M., and A.W. Donaldson, 1950. Use of a surface-active agent to prevent transfer of malarial parasites between blood films during mass staining procedures. *J. Parasitol.* 36:84.
3. Melvin, D. M., and M.M. Brooke, 1955. Triton X-100™ in Giemsa staining of blood parasites. *Stain Technol.* 30:269-275.
4. U.S. Dept. of Health and Human Services, Centers for Disease Control and Prevention.
5. Wilcox, A., 1960. *Manual for the Microscopical Diagnosis of Malaria in Man*. U.S. Public Health Service publication no. 796. U.S. Government Printing Office, Washington, D.C. (Descontinuado)

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**Giemsa (Malaria) Stain
Methylene Blue Phosphate Pre-Stain
Giemsa (Malaria) Stain Buffer**

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Manufactured by CalibreScientific US, Inc.
1311 SE Cardinal Court, Suite 170
Vancouver, WA 98683 USA

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QC1™ est une marque de commerce d' CalibreScientific AMER IP LLC., 6201 Trust Dr, Holland, OH 43528
MDSS GmbH