Auramine O/Rhodamine B Fluorescent Stain Set, Auramine O/Rhodamine B Fluorescent Stain, TB Fluorescent Decolorizer, Potassium Permanganate Counterstain
TABLE OF CONTENTS

English...	1	German...	7
French...	3	Product Codes...	9
Italian...	5	Glossary of Symbols...	9

INTENDED USE

Alpha-Tec Systems, Inc. Auramine O Fluorescent Stain Set and components are recommended for use in qualitative procedures for the fluorescent microscopic detection of mycobacteria.

SUMMARY

One of the earliest methods devised for the detection of the tubercle bacillus is the microscopic staining technique. Mycolic acid within the mycobacteria cell wall interacts with dyes resulting in the characteristics known as "acid-fastness". Acid-fast microscopy is the most rapid, initial step in diagnosis and provides information about the number of acid-fast bacilli present. Hagermann described the use of fluorescent dyes for the detection of acid-fast bacilli in clinical specimens. In 1962, Truant, Brett, and Thomas evaluated the usefulness of a fluorescent staining technique and found it to yield a larger number of positive smears than the conventional fuchsin-stained method. In 1966, Bennedson and Larsen also reported a higher yield of positive smears with a reduced time requirement when using the fluorescent technique.

The fluorochrome dyes, Auramine O/Rhodamine B, used in this stain reacts with the mycolic acids in the acid-fast cell wall of the organism and is refractory when rinsed with acid-alcohol (TB Fluorescent Decolorizer). The counterstain, Potassium Permanganate, renders tissue and debris non-fluorescent thereby reducing the possibility of artifacts. The cells visualized under ultraviolet light appear bright reddish-orange.

FOR IN VITRO DIAGNOSTIC USE ONLY
PRECAUTIONS

Auramine O/Rhodamine B Fluorescent Stain is poisonous and may be harmful or fatal if swallowed. May cause irritation to skin, eyes and respiratory tract. Combustible, keep away from heat and flame. Contains material which may cause cancer based on animal data. Avoid breathing vapor and eye/skin contact. TB Fluorescent Decolorizer may be harmful if swallowed. Avoid breathing vapor and eye/skin contact. Flammable. Keep away from heat sources.

STORAGE

This product is ready for use, and no further preparation is necessary. Store in its original container at 15-30°C until used. Store away from heat and direct sunlight.

USER QUALITY CONTROL

This product should not be used if the color has changed, the expiration date has passed, or there are other signs of deterioration.

SPECIMEN COLLECTION AND PREPARATION

Appropriate specimens for the detection of *Mycobacterium* spp. should be collected according to prescribed standards and delivered to the laboratory in a safe and timely manner. Refer to local procedural guidelines for this information.

PROCEDURE

Materials Provided: Auramine O/Rhodamine B Fluorescent Stain, TB Fluorescent Decolorizer, Potassium Permanganate Counterstain. Each component is also available individually.

Materials Not Provided: Pipettes, glass slides, fluorescent microscope, immersion oil, deionized or distilled water, Bunsen burner or slide warmer, QC1™ AFB Stain Control Slide (#0003240).

SPECIMEN PROCESSING

1. Make a thin smear of the specimen over an area of 2-3 sq cm. Heat fix by warming on a slide warmer (65-75°C) for 10 minutes, or until the smear is dry. If a Bunsen burner is used, pass the slide slowly 3 times over its cone of heat. Do not scorch.
2. Flood the smear with Auramine O Fluorescent Stain for 15 minutes.
3. Rinse with deionized or distilled water and drain.
4. Decolorize with TB Fluorescent Decolorizer for 2-3 minutes.
5. Rinse with deionized or distilled water and drain.
6. Flood smear with Potassium Permanganate Counterstain for 2-4 minutes, no longer.
7. Rinse with deionized or distilled water and allow to air-dry.
8. Examine microscopically using low power (10x – 40x) under a fluorescent microscope and confirm at 100x under oil immersion.
9. A positive fluorescent smear may be restained by the conventional Ziehl-Neelsen or Kinyoun procedure.

EXPECTED RESULTS

Positive Reaction: Acid-fast positive organisms fluoresce bright reddish-orange against a dark background.

Negative Reaction: Nonacid-fast organisms will not fluoresce or may appear a pale red, quite distinct from the bright acid-fast organisms.

LIMITATIONS OF PROCEDURES

1. Overheating (burning) during fixation should be avoided.
2. Do not stain smears which have only been air-dried. Smears must also be "fixed".
3. Smears should not be too thick. After air-drying, examine under a microscope.
4. After staining it is essential that the back surface of the slide is wiped clean.
5. If washing with deionized or distilled water is not done adequately, crystallization of the stain may appear on the slide.
6. A positive staining reaction provides presumptive evidence of the presence of mycobacteria. A negative staining reaction does not indicate that the specimen will be culturally negative. Therefore, cultural methods must be employed.
7. Most strains of rapid growers may not appear fluorescent. It is recommended that all negative fluorescent smears be confirmed with Kinyoun Carbofuchsin Stain (#0003524), Ziehl-Neelsen (#0003531), or Rapid TB Stain (#0003445). At least 100 fields should be examined before being reported as negative.
8. Excessive exposure to the counterstain may result in a loss of brilliance of the fluorescing organism.
9. Turbidity may develop in the stain, but it will not interfere with its effectiveness. Shake the bottle before using.
10. Stained smears should be observed within 24 hours of staining because of the possibility of the fluorescence fading.

SPECIFIC PERFORMANCE CHARACTERISTICS

This product has been tested using the QC1™ AFB Quality Control Slide and has been found to yield acceptable stain results. Appropriate testing should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

BIBLIOGRAPHY

1. Baron, E.J., L.R. Peterson, and S.M. Finegold. 1994. Bailey and Scott's Diagnostic Microbiology. 9th Ed. Mosby, St. Louis, MO.
2. Bennedson, J. and S.O. Larsen. 1966. Scand, J. Respir. Dis. 47:114-120.
3. Hagermann, P.K.H. 1937. Zbl. Bakt. 140:184.
4. Kent, P.T. and G.P. Kubica. 1985. Public Health Mycobacteriology. A Guide for the Level III Laboratory. U.S. Dept. of Health and Human Services. Centers for Disease Control, Atlanta, GA.
5. Kinyoun, J.J., 1915, Am. J. Public Health. 5:867-870.
6. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover. 1995. Manual of Clinical Microbiology. 6th Ed. ASM, Washington, D.C.

[+1] 800.221.6058 United States
[+1] 360.260.2779 International

Directions For Use for the following:

Auramine O/Rhodamine B Fluorescent Stain Set, Auramine O/Rhodamine B Fluorescent Stain, TB Fluorescent Decolorizer, Potassium Permanganate Counterstain

7. Truant, J.P., W.A. Brett, and W. Thomas. 1962. Henry Ford Hosp. Med. Bull. 10:287-296.

CONTACT

Alpha-Tec Systems, Inc. offers a complete line of reagents, stains, and QC1™ Quality Control Slides for AFB, Parasitology, Bacteriology, and Mycology processing, as well as O&P collection systems and concentration devices for Parasitology. For Technical Assistance, email Technical@AlphaTecSystems.com, and for Customer Service, email Sales@AlphaTecSystems.com, or call either [+1] 800.221.6058 (USA) or [+1] 360.260.2779 between 8AM and 4PM Monday through Friday, Pacific Time.

WARRANTY

This product is warranted by Alpha-Tec Systems, Inc. to perform as described in the labeling and literature supplied. Alpha-Tec Systems, Inc. disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall Alpha-Tec Systems, Inc. be liable for any consequential damages arising out of aforesaid express warranty.

TRADEMARKS

QC1™ is a trademark of Alpha-Tec Systems, Inc., 1311 SE Cardinal Court, Suite 170, Vancouver, WA 98683 USA.

Auramine O/Rhodamine B Fluorescent Stain Set, Auramine O/Rhodamine B Fluorescent Stain, TB Fluorescent Decolorizer, Potassium Permanganate Counterstain

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QC1™ ist ein Warenzeichen von Alpha-Tec Systems, Inc., 1311 SE Cardinal Court, Suite 170, Vancouver, WA 98683 USA.
Manufactured by Alpha-Tec Systems, Inc.
MDSS GmbH