TABLE OF CONTENTS

English...	1	Product Codes...	13
Spanish...	5	Glossary of Symbols...	13
Italian...	9		

INTENDED USE

PROTO-FIX® fixative is designed for the collection, transport, and preservation of fecal specimens for permanent staining and concentration procedures to diagnose intestinal parasites.

PROTO-FIX fixative is a clear, colorless, one-vial processing fixative for wet preparations, permanent stains, concentrations, DFA, EIA, and PCR methodologies for fecal specimens. PROTO-FIX is a low-alcohol, low-viscosity fixative containing no heavy metals, PVA, or aniline dyes.

SUMMARY

Intestinal parasitic infections are diagnosed by recovery and identification of protozoan trophozoites and/or cysts or helminth eggs and/or juvenile nematodes. Treatment is based upon positive identification of pathogenic species. As circumstances often prevent immediate examination of fecal specimens, the specimen must be preserved for later examination.

PROTO-FIX fixative is an alternative to the traditional two-vial system of PVA metal fixatives and formalin. PROTO-FIX fixed smears can be stained using modified Wheatley's Trichrome Stain. Protozoan trophozoites, cysts, helminth eggs, and juvenile nematodes can be concentrated using the CONSED® Concentration Procedure of the FEA [formalin ethyl acetate] procedure. The CONSED Concentration Reagent in conjunction with the performance characteristics of PROTO-FIX yield equal or better concentration and staining results than conventional methodologies (PVA and 10% formalin) and provides the convenience and cost savings of a single-vial system. PROTO-FIX fixative has no PVA, mercury, zinc or other heavy metals.

The CONSED Concentration Procedure concentrates all stages of all parasites present in the fecal specimen, including protozoan trophozoites. Trophozoites are clearly visible and identifiable on slides prepared from the concentrated specimen. PROTO-FIX, when used with CONSED, will also concentrate higher numbers of juvenile nematodes [*Strongyloides* spp.]. PROTO-FIX fixative, unlike 10% formalin and SAF, kills fertilized eggs of *Ascaris* spp. ELISA procedures and DFA procedures can be performed directly from fecal specimens collected and fixed in PROTO-FIX.

FOR IN VITRO DIAGNOSTIC USE ONLY

PRECAUTIONS

1. PROTO-FIX is poisonous. If ingested, provide plenty of milk or water. Contact a physician or poison center immediately.
2. Avoid contact with skin and eyes. If contacted, wash thoroughly with water. Contact a physician if irritation occurs.
3. Observe all safety precautions for handling stool specimens.

STABILITY AND STORAGE

Store at room temperature (15-30°C). Avoid excessive heat and sunlight.

USER QUALITY CONTROL

PROTO-FIX fixative is a colorless, clear precipitate-free liquid. If excessive precipitation occurs, discard the fixative vial. The fixative should not be used beyond the product expiration date. The use of a positive parasite control slide for permanent staining procedures is recommended. The PROTO-FIX control slides (#0003255) contain *G. lamblia* trophozoites, ready for staining.

SPECIMEN COLLECTION AND PREPARATION

1. The patient should be instructed not to take antacids, oily laxatives or antidiarrheal medication unless prescribed by a physician prior to collection of the sample(s). Radiological examination utilizing

contrast chemicals (i.e. bismuth or barium) should be avoided prior to the collection of the fecal specimen for parasitic analysis, as they block the staining of ova and parasites.

2. Multiple fecal specimens are recommended to increase the predictive value of parasite recovery.
3. Fecal specimens should be passed into a dry, clean container. Do not contaminate the feces with water or urine.
4. Add enough specimen to each vial to reach the 'FILL LINE' found on the label of the vial. When present, slimy, watery or bloody portions should be collected.
5. Mix specimens well. Replace the cap and shake the vial vigorously.
6. Label the vials with appropriate patient information and return them to the laboratory. Do not freeze or incubate the specimen/fixative vials. Keep the specimen/fixative from heat.
7. After mixing thoroughly, specimen must be fixed for a minimum of one hour to assure adequate fixation of the parasitic elements.

PROCEDURE

Material Provided:

PROTO-FIX is available in various fill volumes. PROTO-FIX may be packaged alone, in a single patient vial set, or in a two-vial patient collection set (clean vial or ETM®).

Materials Not Provided:

Pasteur pipettes, centrifuge with a free-swinging head, microscope, microscope slides and coverslips, cotton-tipped applicators, reagent, and materials for staining and concentration are not provided.

Recommended Additional Reagents and Materials:

CELL-BOND® microscope slides (#0003257), CONSED (#0004628), Ethyl Acetate (#0003344), PRS™ (0004044) or PARA-PRO® *fc50* Concentration System (#0004060) are recommended.

MACROSCOPIC EXAMINATION

The fecal specimen should be examined for consistency. Record any visible irregularities such as worms, proglottids, mucus and/or blood.

DIRECT MICROSCOPIC EXAMINATION

The use of microscopic examination holds limited diagnostic value and should be limited to prescreening of specimens in the field when a centrifuge is not available. If the laboratory is going to perform the concentration procedure and/or the permanent stain, the direct microscopic examination can be eliminated.

CONCENTRATION PROCEDURES

NOTE: The CONSED Concentration Reagent is recommended as it increases the recovery of ova, helminths, and parasites. In addition, the CONSED Concentration Procedure increases the laboratory efficiencies and the diagnostic value of the permanent stain, as the permanent smear can be performed from the CONSED concentrated pellet.

1. Thoroughly mix the PROTO-FIX fixed specimen by shaking the specimen/fixative vial.
2. Add 4-5 drops of Triton X-100 to the specimen/fixative vial (up to 8 drops may be added if the specimen is highly mucoid).
3. Re-cap the fixative vial and mix the contents thoroughly by shaking for 10-20 seconds.
4. Place a PRS concentration funnel or a PARA-PRO *fc50* concentration unit (or other filtration device) into an appropriate polypropylene centrifuge tube. Remix the fixed specimen and pour the sample through the funnel into the receiver centrifuge tube.

NOTE: do not force the fixed specimen solution through the funnel system.

5. Following the filtration process, retain the receiver centrifuge tube and discard the filtration unit into an appropriate disposal container.
6. Transfer 2 ml of the filtered specimen into a 15 ml centrifuge tube. To the 2 ml of the filtered PROTO-FIX fixed specimen, add 8 ml of the CONSED reagent and 4 ml of ethyl acetate (or replacement reagent) to the sample in the centrifuge tube. Replace the cap on the tube, invert the tube, and shake vigorously for 30 seconds.

Pressure will build up within the tube during shaking. To remove this pressure, loosen the cap carefully, and then retighten the cap prior to centrifugation.

7. Place the capped centrifuge tubes into the centrifuge (with a free-swinging head) and centrifuge for 10 minutes at 500 to 600 xg. Following centrifugation, four layers will develop:
 - (1) A top layer of mostly ethyl acetate
 - (2) An interface layer of fatty fecal debris
 - (3) A lower solution layer
 - (4) A pellet/sediment layer
8. Holding the centrifuge tube in a vertical position, remove the cap and free the plug of debris from the sides of the tube by ringing the tube with a wooden applicator stick. Carefully pour the top three layers into an appropriate waste container. **NOTE:** if the pellet begins to break up, place the tube in an upright position, save the pellet, and carefully aspirate any residual reagent off of the pellet with a pipette). While the tube is still tipped in the decanting position, use cotton-tipped swabs to remove remaining debris and ethyl acetate from the sides of the tube. Do not turn the tube upright until the sides of the tube have been thoroughly cleaned of the reagent solutions. If the sides are not cleaned thoroughly with the swab, lipid droplets can mix with the sediment pellet making the microscopic examination more difficult. Allowing excess ethyl acetate to run back into the pelleted sediment will result in a poor wet mount preparation due to the formation of solvent bubbles.
9. Add 3-6 drops (or an amount equal to the volume of the pellet) of PROTO-FIX. Using an applicator stick, mix thoroughly. (**NOTE:** The smear for the permanent stain and slides for special stains can be made at this point in the procedure. See the section header "Preparing Slides for Smears", and "Miscellaneous Procedures").
10. Prepare the wet mount by placing 1 drop of the pellet prepared in step #9 onto a clean glass slide. Add a drop of iodine solution (Lugol's Iodine or Dobell & O'Connor's Iodine) to the specimen, mix gently, and coverslip. Examine the slide microscopically for ova, helminths, and parasites. Consult appropriate references for the identification of ova and parasites.

FORMALIN/ETHYL ACETATE CONCENTRATION PROCEDURE

1. Thoroughly mix the PROTO-FIX fixed specimen by shaking the specimen/fixative vial.
2. Add 4-5 drops of Triton X-100 to the specimen/fixative vial. (Up to 8 drops may be added if the specimen is highly mucoid.)
3. Re-cap the fixative vial and mix the contents thoroughly by shaking for 10 to 20 seconds.
4. Place a PRS concentration funnel or a PARA-PRO *fc50* concentration unit (or other filtration device) into an appropriate polypropylene centrifuge tube. Remix the fixed specimen and pour the sample through the funnel into the receiver centrifuge tube. **NOTE:** Do not force the fixed specimen solution through the funnel system.
5. Following the filtration process, retain the receiver centrifuge tube and discard the filtration unit into an appropriate disposal container.
6. Transfer 3 ml of the filtered specimen into a 15 ml centrifuge tube. To the 3 ml of the filtered PROTO-FIX fixed specimen, add 7 ml of 10% buffered formalin and mix the specimen.
7. Add 4 ml of ethyl acetate (or replacement reagent) to the sample in the centrifuge tube. Replace the cap on the tube, invert the tube, and shake vigorously for 30 seconds. Pressure will build up within the tube during shaking. To remove this pressure, loosen the cap carefully and retighten the cap prior to centrifugation.
8. Place the capped centrifuge tubes into the centrifuge (with a free swinging head) and centrifuge for 10 minutes at 500 to 600 xg. Following centrifugation four layers will develop:
 - (1) A top layer of mostly ethyl acetate
 - (2) An interface layer of fatty fecal debris
 - (3) A lower solution layer
 - (4) A pellet/sediment layer
9. Holding the centrifuge tube in a vertical position, remove the cap, free the plug of debris from the sides of the tube by ringing the tube with a wooden applicator stick. Carefully pour the top three layers

into an appropriate waste container. **NOTE:** If the pellet begins to break up, quickly upright the tube to save the pellet and carefully aspirate any residual reagent off of the pellet with a pipette). While the tube is still tipped in the decanting position, use cotton-tipped swabs to remove remaining debris and ethyl acetate from the sides of the tube. Do not turn the tube upright until the sides of the tube have been thoroughly cleaned of the reagent solutions. If the sides are not cleaned thoroughly with the swab, lipid droplets can mix with the sediment pellet making the microscopic examination more difficult. Allowing excess ethyl acetate to run back into the pelleted sediment will result in a poor wet mount preparation due to the formation of solvent bubbles.

10. Add a few drops of PROTO-FIX to the pellet and mix well.
11. Prepare the wet mount by placing 1 drop of the sediment prepared in step #10 onto a clean glass slide. Add a drop of iodine solution (Lugol's Iodine or Dobell & O'Connor's Iodine) to the specimen, mix gently, and coverslip. Examine the slide microscopically for ova, helminths, and parasites. Consult appropriate references for the identification of ova and parasites.

PREPARING SLIDES FOR SMEARS

PROTO-FIX contains no PVA and therefore requires the use of an adhesive to improve the adhesion of the specimen to the slide during routine staining procedures. Use the following method in the slide preparation:

CELL-BOND® Slides are ready-to-use bio adhesive impregnated slides, which contain a charge, binding the specimen to the slide. Unlike slide Coating Solutions, Mayer's Albumin or PVA, CELL-BOND Slides provide not competitive staining background. This results in a clearer microscopic view of the fecal smear.

PERMANENT SMEAR PREPARATION

1. Transfer 1-2 drops of the CONSED concentrated pellet or 1-2 drops of the unconcentrated PROTO-FIX fixed specimen to a CELL-BOND slide (see above).
2. Gently and evenly spread the sample over the microscope slide. Lay or hold the slide flat with the specimen side up. Using an applicator stick, chop the specimen, spreading specimen out. This will create thick and thin areas.
3. Allow the slide to remain flat for 1 to 2 minutes. If there is still excessive liquid on the slide, stand the slide in a drying rack at 45 degrees to allow any excess liquid to drain. **NOTE:** Slides should dry for 10-15 minutes. Fecal films/smears dry at different rates. If some fecal smears are still wet after 15 minutes, they may be air-dried with a cool fan or blower. (Do not use heat to dry the fecal smear slides). Fecal films/smears can be slightly wet when staining begins. When the excess liquid stops draining (usually about 5 to 10 minutes, carefully wipe away any excess liquid on the end or edges of the slide and place the slides into a carrier for staining.

PERMANENT STAIN PROCEDURE

Modified Wheatley's Trichrome Stain Procedure

Reagent	Timing
70% Ethanol (#0003359)	1.5 minutes. Remove and drain off excess liquid by touching edge of slide to absorbent material
70% Ethanol	1.5 minutes. Remove and drain off excess liquid by touching edge of slide to absorbent material
Trichrome Stain (#0003351)* (Wheatley's)	13 minutes. Remove and drain off excess liquid by touching edge of slide to absorbent material.
90% Acid-Ethanol (#0003350)	2 seconds. Remove and drain off excess liquid by touching edge of slide to absorbent material.
100% Ethanol (#0003303)	1 minute. Remove and drain off excess liquid by touching edge of slide to absorbent material.

(Continued on next page)

Reagent	Timing
100% Ethanol (#0003359)	1 minute. Remove and drain off excess liquid by touching edge of slide to absorbent material.
Xylene (#0003342) or PRO-CLEAR™ (#0003336)	3 minutes. Remove and drain off excess liquid by touching edge of slide to absorbent material and coverslip with mounting fluid. Read under oil immersion.

NOTE: Staining time can be varied depending on the intensity desired for the final stain.

NOTE: The xylene substitute AmeriClear does not yield appropriate stain results and should not be used.

MISCELLANEOUS PROCEDURES

- Place 1 drop of the specimen onto a clean glass slide. Using an applicator stick, spread the specimen evenly to achieve a thin layer on the slide. Allow the slide to dry for 5 to 10 minutes before proceeding.
- Dip the dried smear into 70% ethanol for 1 minute.
- Remove the smear from the ethanol and gently shake off any excess ethanol. Dip the smear into distilled water about 8 times to sheet off the ethanol. Shake off the excess water and stain with the Cryptosporidium Stain (#0003357) or other appropriate stains.

DFA Staining Preparation

- Follow the smear preparation procedure listed above under “**Cryptosporidium Staining Procedure**” to prepare a slide for DFA staining.
- Follow the specific DFA manufacturer’s directions for staining the smear for DFA *Cryptosporidium parvum* and *Giardia lamblia*.
- PROTO-FIX has been validated against the MERIFLUOR DFA stain for *Cryptosporidium parvum* and *Giardia lamblia*.

Microsporidium Staining Procedure

- Follow the smear preparation procedure listed above under “**Cryptosporidium Staining Procedure**” to prepare a slide for Microsporidium staining.
- Follow the specific Microsporidium Stain manufacturer’s directions for staining the smear

Enzyme Immunoassay (EIA) Procedure

PROTO-FIX has been validated against the following manufacturers’ products for EIA testing for *Cryptosporidium parvum* and *Giardia lamblia*.

- ProSpecT *Cryptosporidium* Rapid Assay
- ProSpecT *Cryptosporidium* Microplate Assay
- Wampole CRYPTOSPORIDIUM TEST
- Xpect *Cryptosporidium* Kit
- ProSpecT *Cryptosporidium parvum* and *Giardia lamblia* Microplate ELISA Assay
- COLORPAC *Giardia/Cryptosporidium* Rapid Test
- IMMUNOCARD stat! *Cryptosporidium/Giardia*
- Xpect *Giardia* and *Cryptosporidium* Kit
- ProSpecT *Giardia/Cryptosporidium* Microplate Assay
- Wampole GIARDIA TEST
- ProSpecT *Giardia* Rapid Assay
- Xpect *Giardia* Kit

EXPECTED RESULTS

Using the modified Wheatley’s Trichrome Stain procedure with organisms fixed in PROTO-FIX, the cytoplasm of trophozoites and cysts is blue-green tinged with purple. The nuclear chromatin, chromatid bodies, erythrocytes and bacteria stain purple to purplish-red. Background and artifacts stain green to purple. Helminth eggs and larvae stain red to purplish-red.

LIMITATIONS OF PROCEDURES

- The proper fixation of intestinal ova and parasites yields quality microscopic morphology and staining. Specimens not properly fixed (specimens delayed before being fixed, improper ratios of specimen to fixative, and improper mixing of the specimen into the fixative) may yield poor microscopic morphology making it difficult or impossible to properly identify the ova or parasite. False negative examinations may occur if too little specimen or if too much specimen is used in the concentration procedures.
- The Modified Wheatley’s Trichrome stain referred to in this Directions For Use is designed for use with PROTO-FIX only. All proficiency specimens received for staining (usually Zn-PVA smears) must be stained using the modified Trichrome stain designed for the zinc PVA fixative (#0003125).
- The xylene substitute AmeriClear does not yield appropriate stain results and should not be used.

SPECIFIC PERFORMANCE CHARACTERISTICS

PROTO-FIX fixative, when used according to the procedures listed above, should yield quality microscopic morphologies from the wet preparation permanent stain concentration procedures.

BIBLIOGRAPHY

- Allen, K., Frankel, J.W., and Kwa, B. 1997. “Comparison of CONSED and Formalin-ethyl-Acetate Methods for Concentrating Intestinal Parasites and Eggs.” Southeastern Branch, ASM, 1997 Meeting.
- Amin, O. 2000. “Evaluation of a New System for the Fixation, Concentration and Staining of Intestinal Parasites in Fecal Specimens, with Critical Observation on the Trichrome Stain.” Journal of Microbiological Methods. 38:127-132.
- Burrows, R.B. 1965. “Microscopic Diagnosis of the Parasites of Man.” Yale University Press, New Haven, CT.
- Erdman, D. 1981. “Clinical Comparison of Ethyl Acetate and Diethyl Ether in the Formalin-ether Sedimentation Technique.” J Clin. Micro. 14:483-485.
- Faust, E.C., et al. 1938. “A Critical Study of Clinical Laboratory Techniques for the Diagnosis of Protozoan Cysts and Helminth Eggs in Feces.” Am.J Trop. Hyg. Med. 18:169-183.
- Garcia, L.S., and Laurence Ash. 1979. Diagnostic Parasitology. C.V. Mosby Co., St. Louis, MO, pp. 9-24.
- Garcia, L.S., M. Voge. 1990. Diagnostic Clinical Parasitology. Am. J. Med. Technol. 46:459-467.
- Jensen, B., et al. 2000. “Comparison of Polyvinyl Alcohol Fixative with Three Less Hazardous Fixatives for the Detection and Identification of Intestinal Parasites.” J. Clin. Micro. 38(4):1592-1598.
- Markell, E., and Marietta Voge. 1981. Medical Parasitology. W.B. Saunders Co., Philadelphia, PA, pp. 317-355.
- Melvin, D.M. and M.M. Brooke, 1982. “Laboratory Procedures for the Diagnosis of Intestinal Parasites.” 3rd Edition. US Department of Health and Human Services Pub. No. (CDC) 82-8282, Atlanta, GA.
- Price, D.L. 1994. Procedure Manual for the Diagnosis of Intestinal Parasites. CRC Press, Boca Raton, FL.
- Sholten, T.H. and J. Yang. 1974. “Evaluation of Unpreserved and Preserved Stools for the Detection and Identification of Intestinal Parasites.” Am. J. Clin. Pathol. 62:563-567.
- Williams, R.W. and Dennis, M. 1998. “Technical Performance of a New Fecal Fixative for *Giardia lamblia* and *Cryptosporidium parvum* in Fecal Specimens.” ASM Meeting, Atlanta, GA.
- Young, K., S. Bullock, D. Melvin and C. Spruill. 1979. “Ethyl Acetate as a Substitute for Diethyl Ether in the Formalin-ether Sedimentation Technique.” J. Clin. Micro. 10:852-853.
- Unpublished data on file.

CONTACT

CalibreScientific US, Inc. offers a complete line of reagents, stains, and QC1™ Quality Control Slides for AFB, Parasitology, Bacteriology, and Mycology processing, as well as O&P collection systems and concentration devices for Parasitology. For Technical Assistance, email Technical@AlphaTecSystems.com, and for Customer Service, email Sales@AlphaTecSystems.com, or call either [+1] 800.221.6058 (USA) or [+1] 360.260.2779 between 8AM and 4PM Monday through Friday, Pacific Time.

WARRANTY

This product is warranted by CalibreScientific US, Inc. to perform as described in the labeling and literature supplied. CalibreScientific US, Inc. disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall CalibreScientific US, Inc. be liable for any consequential damages arising out of aforesaid express warranty.

TRADEMARKS

CELL-BOND®, CONSED®, ETM®, PARA-PRO® *fc50*, PRO-CLEAR™, PROTO-FIX®, PRS™, and QC1™ are trademarks of CalibreScientific AMER IP LLC., 6201 Trust Dr, Holland, OH 43528.

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1. Allen, K., Frankel, J.W., and Kwa, B. 1997. "Comparison of CONSED Red and Formalin-ethyl-Acetate Methods for Concentrating Intestinal Parasites and Eggs." Southeastern Branch, ASM, 1997 Meeting.
2. Amin, O. 2000. "Evaluation of a New System for the Fixation, Concentration and Staining of Intestinal Parasites in Fecal Specimens, with Critical Observations on the Trichrome Stain." Journal of Microbiological Methods. 38:127-132.
3. Burrows, R.B. 1965. "Microscopic Diagnosis of the Parasites of Man." Yale University Press, New Haven, CT.
4. Erdman, D. 1981. "Clinical Comparison of Ethyl Acetate and Diethyl Ether in the Formalin-ether Sedimentation Technique." J Clin. Micro. 14:483-485.
5. Faust, E.C., et al. 1938. "A Critical Study of Clinical Laboratory Techniques for the Diagnosis of Protozoan Cysts and Helminth Eggs in Feces." Am. J Trop. Hyg. Med. 18:169-183.
6. Garcia, L.S., and Lawrence Ash. 1979. Diagnostic Parasitology. C.V. Mosby Co., St. Louis, MO, pp. 9-24.
7. Garcia, L.S., M. Voge. 1990. Diagnostic Clinical Parasitology. Am. J. Med. Technol. 46:459-467.
8. Jensen, B., et al. 2000. "Comparison of Polyvinyl Alcohol Fixative with Three Less Hazardous Fixatives for the Detection and Identification of Intestinal Parasites." J. Clin. Micro. 38(4):1592-1598.
9. Markell, E., and Marietta Voge. 1981. Medical Parasitology. W.B. Saunders Co., Philadelphia, PA, pp. 317-355.
10. Melvin, D.M. and M.M. Brooke, 1982. "Laboratory Procedures for the Diagnosis of Intestinal Parasites." 3rd Edition. US Department of Health and Human Services Pub. No. (CDC) 82-8282, Atlanta, GA.
11. Price, D.L. 1994. Procedure Manual for the Diagnosis of Intestinal Parasites. CRC Press, Boca Raton, FL.
12. Sholten, T.H. and J. Yang. 1974. "Evaluation of Unpreserved and Preserved Stools for the Detection and Identification of Intestinal Parasites." Am. J. Clin. Pathol. 62:563-567.
13. Williams, R.W. and Dennis, M. 1998. "Technical Performance of a New Fecal Fixative for *Giardia lamblia* and *Cryptosporidium parvum* in Fecal Specimens." ASM Meeting, Atlanta, GA.
14. Young, K., S. Bullock, D. Melvin and C. Spruill. 1979. "Ethyl Acetate as a Substitute for Diethyl Ether in the Formalin-ether Sedimentation Technique." J. Clin. Micro. 10:852-853.
15. Unpublished data on file.

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MARCAS REGISTRADAS:

CELL-BOND®, CONSED®, ETM®, PARA-PRO® *fc50*, PRO-CLEAR™, PROTO-FIX®, PRS™, and QC1™ son marcas registradas de CalibreScientific AMER IP LLC., 6201 Trust Dr, Holland, OH 43528.

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Manufactured by CalibreScientific US, Inc.
MDSS GmbH