PROTO-FIX®

Trichrome Stain Set (Wheatley's Modification)

TABLE OF CONTENTS

English...	1	Product Codes...	5
Spanish ...	3	Glossary of Symbols...	5

INTENDED USE

The PROTO-FIX® Wheatley's Trichrome Stain Set is a complete set of reagents for use in performing the Wheatley's Trichrome Stain procedure for the identification of intestinal protozoa in fecal specimens fixed in PROTO-FIX.

SUMMARY

Intestinal parasitic infections are diagnosed by recovery and identification of protozoan trophozoites and/or cysts or helminth eggs and/or juveniles. The detection and correct identification of intestinal protozoa is frequently dependent on the examination of a permanently stained smear as smaller protozoa are often missed with only the direct smear and concentration methods. The Trichrome stain has been used since 1929 as a histological stain for muscle tissue. In 1949, Gomori developed a shortened and rapid method for trichrome staining of histologic and cytologic sections. In 1951, Wheatley modified the Gomori procedure and, using Trichrome Stain, developed a rapid staining procedure for intestinal amoeba and flagellates.

FOR IN VITRO DIAGNOSTIC USE ONLY

PRECAUTIONS

This product should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, and media after their use. Directions should be read and followed carefully. Refer to the Safety Data Sheets for additional information.

STORAGE

Store the product in the original container at room temperature (15-30°C). Avoid extremes of temperature and light. Keep all containers tightly closed when not in use.

USER QUALITY CONTROL

A positive parasite smear should be processed with each batch of slides or per the regulatory agency used by your laboratory to verify the quality of the staining reagents and the technique of the procedure. This product should not be used if the color has changed, the expiration date has passed or there are other signs of deterioration.

SPECIMEN COLLECTION AND PREPARATION

Permanent smears can be prepared from the unconcentrated specimen fixed in PROTO-FIX or from a specimen in PROTO-FIX that has been concentrated following the CONSED® Concentration procedure. (Refer to the CONSED directions for use). Specimens should be collected and handled following recommended guidelines.

Preparation of Smears for Staining Specimens Preserved in PROTO-FIX:

PROTO-FIX contains no PVA and as such may require the use of a coated slide, such as a CELL-BOND® Slide (#0003257), to improve the adhesion of the specimen to the slide during staining procedures.

Transfer 1-2 drops of the PROTO-FIX fixed specimen to a CELL-BOND slide). Lay or hold the slide flat with the specimen side up. Using an applicator stick, gently and evenly spread the sample over the slide then, using a chopping motion, spread the specimen out to create thick and thin areas.

Allow the slide to remain flat for 1 to 2 minutes. Fecal smears can be slightly wet when staining begins. If there is still excessive liquid on the slide, stand it in a drying rack at 45° to allow any excess liquid to drain. If fecal smears are still wet after 15 minutes, they may be air-dried with a cool fan (do not use heat to dry the fecal smear slides). When the

excess liquid stops draining, carefully wipe away any excess liquid on the edges of the slide and proceed with staining.

PROCEDURE

The PROTO-FIX Stain Set includes all reagents necessary to perform the modified Wheatley's Trichrome stain procedure with PROTO-FIX specimens. All reagents are packaged in convenient, wide-mouth jars which can be used for the staining procedure. It is not necessary to transfer reagents to Coplin jars.

Materials Provided:

70% Ethanol	(2 x 50 ml)
Wheatley's Trichrome Stain	(1 x 50 ml)
90% Acid Ethanol	(2 x 50 ml) (one replacement included)
100% Ethanol	(4 x 50 ml) (one replacement included)
PRO-CLEAR™	(1 x 50 ml)

Materials Not Provided: Mounting medium, CONSED, coverslips, immersion oil, CELL-BOND, applicator sticks, absorbent paper.

SPECIMEN PROCESSING

1. Remove all jars from packaging box and line up sequentially, starting with "Container 1".
2. When ready to stain slide(s), remove caps from jars and place slide directly into the jar.
3. When finished, close lids tightly and store.

NOTE: Less contamination of reagents will occur if a blotting step is included between each staining step. Allow slides to drain between each solution.

Reagent	Timing
70% Ethanol	1.5 minutes. Remove and drain off excess liquid by touching edge of slide to absorbent material.
70% Ethanol	1.5 minutes. Remove and drain off excess liquid by touching edge of slide to absorbent material.
Trichrome Stain	13 minutes. Remove and drain off excess liquid by touching edge of slide to absorbent material.
90% Acid Ethanol	1-3 seconds. Remove and drain off excess liquid by touching edge of slide to absorbent material.
100% Ethanol*	5-10 seconds. Remove and drain off excess liquid by touching edge of slide to absorbent material.
100% Ethanol*	1 minute. Remove and drain off excess liquid by touching edge of slide to absorbent material.
100% Ethanol*	1 minute. Remove and drain off excess liquid by touching edge of slide to absorbent material.
PRO-CLEAR	3 minutes. Remove and drain off excess liquid by touching edge of slide to absorbent material and apply coverslip.
Mount coverslip and examine under oil immersion lens.	

NOTE: Staining time can be varied depending on the intensity desired for the final stain.

*NOTE: 100% Ethanol-based reagent alcohol

NOTES

- Less contamination of reagents will occur if a blotting step is included between each staining step. Allow slides to drain between each solution, and then blot the edge of the slide only.
- The Trichrome stain may become diluted due to carryover of ethanol.

PROTO-FIX®

Trichrome Stain Set (Wheatley's Modification)

- The stain may be restored by leaving the lid off the jar for several hours to allow evaporation of the ethanol.
- The acid ethanol will continue to decolorize as long as it is in contact with the material. The time stated in the procedure for decolorizing should begin when the slide is placed in the acid ethanol solution and end after the slide is rinsed in 100% ethanol.

All components in the PROTO-FIX Wheatley's Trichrome Stain Set are available in larger sizes. Contact the CalibreScientific US, Inc. Sales Department.

EXPECTED RESULTS

Typical staining reactions using Wheatley's Trichrome Stain Set with specimens adequately fixed in PROTO-FIX and stained are listed below:

The nuclear chromatin, chromatoid bodies, ingested erythrocytes and bacteria are purple to red-violet.

The cytoplasm of trophozoites and cysts stains blue tinted with purple. Background material and artifacts stain blue-green to purple.

Helminth ova and larvae stain red to purple but may be more easily identified in a concentrated wet mount.

Macrophages, leukocytes and yeasts show variable staining from green to blue to purple or red.

Special staining procedures are necessary for the identification of *Cryptosporidium*, *Cyclospora* and *Microsporidia*.

LIMITATIONS OF PROCEDURES

Wheatley's Trichrome Staining method is usually a trouble-free procedure when used as directed. Any problems that may occur will usually be one of the following:

Problem: Poor contrast of the chromatin material.

Cause: Over-decolorizing.

Solution: Decolorizing (acid ethanol) requires only a very brief contact, followed by an immediate dip in ethanol.

Problem: Poor staining of the cytoplasm and the nucleus. Degenerate forms that stain weakly.

Cause: Parasitic and cellular elements have degenerated because of improper fixation.

Solution: To ensure proper fixation, specimens must be placed in a fixative solution immediately after passage. The proportion of specimen to fixative (1:3) must be observed, be thoroughly mixed and have sufficient time for fixation.

Problem: Stained preparation is "cloudy" with poor contrast of cellular detail.

Cause: Carryover of solutions from one step to another. Excessive staining (more than 30-40 slides) will weaken or dilute the stain.

Solution: Change all solutions regularly to avoid staining clarity problems

Problem: Inadequate material for microscopic examination.

Cause: Smear is too thin.

Solution: Smears should be moderately thick. If necessary, concentrate the fixed specimen by centrifugation. Pour off excess specimen fixative to provide thicker material for smears.

Problem: Persistent staining failure after considering all of the above.

Cause: Faulty technique and/or contaminated reagents.

Solution: Discard the entire stain series. Try again with new staining solutions.

BIBLIOGRAPHY

1. Balows, A., W.J. Hausler, Jr., K.L. Hermann, H. Isenberg and J.J. Shadomy, Manual of Clinical Microbiology, 5th Edition, ASM, Washington, D.C. 1991.
2. Baron, E.J. and S.M. Finegold, Baily and Scott's Diagnostic Parasitology, 8th Edition. The C.V. Mosby Co. St. Louis, MO 786-790. 1990.
3. Garcia, L.S. and D.A. Bruckner, Diagnostic Medical Parasitology, Elsevier, New York, NY 360-371, 450. 1988.
4. Gormori, G., Am. J. Clin. Pathol., 20:661-663. 1950.
5. Melvin, D.M. and M.M. Brooke, Laboratory Procedures for the Diagnosis of Intestinal Parasites, 3rd Edition, U.S. Dept. of Health and Human Services, Centers for Disease Control, Pub. No. (CDC) 82-8282, 125-129. 1982.
6. National Committee for Clinical Laboratory Standards. Procedures for the Recovery and Identification of Parasites from the Intestinal Tract; Proposed Guidelines. NCCLS Publication M28-P, Villanova, PA 1993.
7. Wheatley, W.B., Am. J. Clin. Pathol., 21:990-991. 1951.

CONTACT

CalibreScientific US, Inc. offers a complete line of reagents, stains, and QC1™ Quality Control Slides for AFB, Parasitology, Bacteriology, and Mycology processing, as well as O&P collection systems and concentration devices for Parasitology. For Technical Assistance, email Technical@AlphaTecSystems.com, and for Customer Service, email Sales@AlphaTecSystems.com, or call either [+1] 800.221.6058 (USA) or [+1] 360.260.2779 between 8AM and 4PM Monday through Friday, Pacific Time.

WARRANTY

This product is warranted by CalibreScientific US, Inc. to perform as described in the labeling and literature supplied. CalibreScientific US, Inc. disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall CalibreScientific US, Inc. be liable for any consequential damages arising out of aforesaid express warranty.

TRADEMARKS:

CELL-BOND®, CONSED®, PRO-CLEAR™, PROTO-FIX®, and QC1™ are trademarks of CalibreScientific AMER IP LLC., 6201 Trust Dr, Holland, OH 43528.

PROTO-FIX®

Trichrome Stain Set (Wheatley's Modification)

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Manufactured by CalibreScientific US, Inc.
MDSS GmbH