

PRINCIPLE

NAC-PAC® RED is used in the qualitative procedure in the N-acetyl-L-cysteine (NALC) digestion and decontamination procedure of clinical specimens for the increased recovery of *Mycobacterium* spp., incorporating phenol red as a pH indicator throughout the decontamination and buffering procedure.

CLINICAL SIGNIFICANCE

The decontamination and digestion procedure, utilizing the compound N-acetyl-L-cysteine (NALC) combined with sodium hydroxide and sodium citrate (trisodium citrate) solution, results in increased yields of tubercle bacilli. The NALC procedure utilizes N-acetyl-L-cysteine as a mucolytic compound by disrupting chemical bonds in mucus. The sodium hydroxide acts as a bacterial decontaminant and the sodium citrate (trisodium citrate) solution stabilizes the NALC by chelating (binding) any heavy metal ions present in the specimen. Since the sodium hydroxide has a pH of approximately 13.00, it will kill bacteria (including mycobacteria after 15-20 minutes of exposure). As such, timing of the decontamination is critical to limit the amount of *Mycobacterium* spp. killed by the basic pH. A pH indicator is incorporated into the decontamination reagent to monitor the pH throughout the decontamination and buffering procedure, allowing the laboratory technologist to visually see when neutralization has been achieved. Bringing the pH to a neutral range can stop the decontamination process. The NPC-67® Neutralizing Buffer or XPR-PLUS® Neutralizing Buffer is used to neutralize the NaOH following the appropriate digestion and decontamination time, resulting in a pH below 8.10. Adding conventional M/15 phosphate buffer or phosphate buffered saline will result in a pH range of 9.40 to 12.20, requiring a titration to a neutral pH with 1N HCL, or continued decontamination of *Mycobacterium* spp. will occur. Studies have documented that pH values above 8.10 are toxic to *Mycobacterium* spp., including *Mycobacterium tuberculosis*. Following the decanting step, PRB™ Pellet Resuspension Buffer is added to achieve a tight neutral pH value (6.80-7.10) in the specimen sediment, optimizing mycobacteria recovery.

SPECIMEN COLLECTION AND PREPARATION

Appropriate specimens for the detection of *Mycobacterium* spp. should be collected according to prescribed standards and delivered to the laboratory in a safe and timely manner. Refer to local procedural guidelines for this information. **FOR IN VITRO DIAGNOSTIC USE ONLY.**

REAGENTS AND MATERIALS

1. **Provided**
 - a. NAC-PAC RED and NALC
2. **Not Provided**
 - a. NPC-67 Neutralizing Buffer
 - b. XPR-PLUS Neutralizing Buffer (for neutralization of digestion/decontamination solutions with \geq 3% NaOH)
 - c. PRB Pellet Resuspension Buffer
 - d. Centrifuge
 - e. Vortex mixer
 - f. Sterile pipettes
 - g. Microscope slides
 - h. TB media
 - i. Centrifuge tubes
 - j. CELL-BOND® Slides.
3. **Storage:** Prior to opening, store at room temperature (15-30° C). After opening, store between 2-8° C. Do not freeze or heat above 30° C. Allow the product to come to room temperature prior to use.
4. **Stability:** NAC-PAC RED is stable to the stated expiration date when stored at the required temperature.

CALIBRATION

N/A

QUALITY CONTROL

Any product showing cloudiness, turbidity, precipitation or coloration should be discarded. Quality controlled microorganisms should be utilized to verify procedures, media and reagents as appropriate for your laboratory's applicable regulatory agency or local procedural guidelines.

PROCEDURE

PRECAUTIONS

The decontamination reagent contains a caustic chemical (sodium hydroxide). Use appropriate care in the handling of this reagent. All clinical specimens submitted for the diagnosis of tuberculosis and other *Mycobacterium* spp. must be treated with appropriate care so as not to contaminate other specimens or laboratory personnel. Use all approved and regulated equipment for processing and detection procedures.

1. Line up specimens (in centrifuge tubes) in a biosafety hood.
2. Loosen specimen container caps. Work in sets equivalent to a centrifuge load.
3. Open the bottle labeled "NAC-PAC RED". Add the NALC powder to the NAC-PAC RED bottle. Shake well to dissolve the NALC powder. **NOTE:** Some residual NALC powder may remain in the vial. It is not necessary to liquefy the portion remaining in the vial. **THIS SOLUTION WILL BE GOOD FOR ONLY 72 HOURS AFTER MIXED.** Discard the mixed solution after 72 hours.
4. To the sterile 50 ml centrifuge tube containing the specimen to be digested, add the NAC-PAC RED / NALC solution in the following amounts:
 - a. For specimens **1-5 ml add a volume of NAC-PAC RED / NALC equal to that of the specimen volume.**
 - b. For specimens **6-7 ml add 5 ml of NAC-PAC RED / NALC.**
 - c. For specimens **8-10 ml add an equal volume of NAC-PAC RED / NALC and split the specimen after step 6** equally into two centrifuge tubes, proceed with steps 7-9 and then combine the sediments from both tubes into one centrifuge tube and proceed with step #10.

Following this protocol will help achieve effective decontamination while also allowing for proper neutralization. If you routinely encounter specimens greater than 10 ml in volume, please contact Alpha-Tec Systems Technical Services for special instructions.

5. Tighten the caps on the centrifuge tubes. Mix each specimen on a vortex until liquefied (30 seconds per specimen).
6. Allow each specimen to stand for 15-20 minutes. Vortex every five minutes during this step.
7. Fill each tube with NPC-67 or XPR-PLUS until effective base pH neutralization is indicated by a color change from red / pink to colorless. Once a colorless point has been reached, do not continue to add NPC-67 or XPR-PLUS to the sample. Tighten cap and swirl by hand to mix. **NOTE:** NPC-67 will achieve a neutral pH (colorless solution) when added to NAC-PAC RED with a NaOH concentration of 3% or lower. For NAC-PAC RED containing 4% NaOH, use XPR-PLUS.
8. Centrifuge the specimen tubes at 3000 xg for 15 minutes. It is recommended but not required to use a refrigerated centrifuge. Each laboratory must check the centrifuge head radius, and use an appropriate nomogram for proper speed selection [rpm] to achieve the desired relative centrifugal field of 3000 xg.
9. Working in a biosafety hood, pour off **all** supernatant into a splash-proof container holding an appropriate disinfectant. Use an appropriate disinfectant to disinfect any contamination on the lip of the specimen tube. Do not allow the disinfectant to run down inside the specimen tube.
10. Resuspend the pellet with 0.5 ml-1.0 ml of PRB. Do not resuspend the pellet with NPC-67, XPR-PLUS, water or saline. **NOTE:** To maximize time to detection for rapid growth automated detection systems, resuspend the pellet with 1.0 ml of PRB. Depending on the needs of your laboratory, the pellet may be resuspended in 0.5 ml of PRB to create a more concentrated sample for increased acid-fast smear sensitivity. Once the smears have been made, add an additional 1.0 ml of PRB to inoculate rapid broth detection systems and other media.
11. Mix the sediment and buffer well, and inoculate the liquid broth for your automated detection equipment per the manufacturer's instructions.
12. Place two drops of the sediment onto the surface of each of the TB media used. **NOTE:** A contamination control plate (BAP or TSA) can be inoculated at this point and incubated at 35-37° C for 48 hours.
13. Make smears for acid-fast staining. Use adhesive CELL-BOND Slides or appropriate sterile albumin adhesive solutions to attach the specimen to the slide. Dry the smears and proceed with acid-fast staining per the manufacturer's directions. **NOTE:** An acid-fast stain control slide should be stained in conjunction with the patient smears to verify the staining technique and components. Call Alpha-Tec Systems, Inc. for a complete list of acid-fast stains and control slides.
14. To the unused portion of the specimen, add the balance of the PRB and refrigerate at 2-8° C to save for future diagnostic procedures or reprocessing if necessary.

CALCULATIONS

N/A

RESULTS

If *Mycobacterium* spp. are present in the clinical specimen and processed according to the procedures listed within this document, the recovery of cultivable, viable, and clinically significant *Mycobacterium* spp. can be expected. NAC-PAC RED was tested on clinical samples and recovered all culture appropriate *Mycobacterium* spp. when the designated procedures were followed.

LIMITATIONS

Timing of the decontamination step, proper buffering, speed and timing of the centrifugation step, proper decanting and addition of the PRB to the pellet are vital to the recovery of *Mycobacterium* spp. Failure to follow the listed procedures may result in decreased numbers of *Mycobacterium* spp. or total loss of *Mycobacterium* spp. resulting in an inaccurate culture report.

NOTES

©2022 Alpha-Tec Systems Inc. All Rights Reserved. CC22-178, Effective Date: 15 Aug 2022, CLSI.L004302.E

Page 2 of 4

1. Procedure Notes

- a. Molecular Diagnostics
NAC-PAC RED has been validated for use with multiple molecular diagnostic methods and systems. For more information regarding compatibility with specific methods or systems, contact Alpha-Tec Technical Services.
- b. Small Volume Specimens
Small volume specimens with correspondingly low post neutralization volumes can make centrifuge balancing difficult. If your laboratory frequently encounters small volume specimens, it is acceptable to add **sterile** saline to the sample to reach a combined volume of 5 ml prior to the addition of NAC-PAC RED / NALC solution. In this case, the sample should be decontaminated with 5 ml of NAC-PAC RED / NALC solution. This will increase the final post neutralization specimen volume making centrifuge balancing easier.
- c. Specimens contaminated with *Pseudomonas* spp.
Specimens contaminated with *Pseudomonas* spp. will need additional treatment with 5% Oxalic Acid (OxA® Oxalic Acid Reagent Kit #0004805). Refer to the Oxalic Acid Directions For Use for complete instructions, or call Alpha-Tec Systems, Inc. Technical Services for information on the pH effects of the Oxalic Acid procedure and the appropriate buffering requirements.
- d. Bloody Specimens
Following the decontamination of the specimen with NAC-PAC RED, bloody specimens may remain pink after the addition of the NPC-67 or XPR-PLUS due to the residual hemoglobin in the specimen. If the color change cannot be visualized due to hemoglobin, add the NPC-67 or XPR-PLUS up to the 50 ml mark to ensure complete neutralization. For additional information, contact Alpha-Tec Technical Services.

2. Summary of Technology

- a. Decontamination and Digestion
 - i. Sodium Hydroxide
 1. Digests bacteria (including *Mycobacterium* spp.) utilizing a high, basic pH.
 2. Mucolytic compound that disrupts chemical bonds in mucus resulting in total specimen digestion.
 - ii. N-acetyl-L-cysteine (NALC)
 1. Mucolytic compound that disrupts chemical bonds in mucus resulting in total specimen digestion.
 2. Combines with sodium hydroxide and trisodium citrate resulting in increased yields of tubercle bacilli.
 - iii. Timing
 1. Timing is critical so as not to limit the die-off of *Mycobacterium* species present in the patient specimen by the basic pH.
- b. Neutralization
 - i. NPC-67 Neutralizing Buffer or XPR-PLUS Neutralizing Buffer
 1. Used to neutralize the NALC reagents following the appropriate digestion decontamination time, resulting in a pH ≤ 8.10 .
 2. Following the decanting step, PRB is added to achieve a tight neutral pH value (6.8-7.1) in the specimen sediment, optimizing *Mycobacteria* recovery.
 3. Studies have documented that pH values above 8.1 are toxic to *Mycobacterium* spp., including *Mycobacterium tuberculosis*. Bringing the pH to a neutral range can stop the digestion procedure.
 4. A pH indicator is incorporated in the digestion decontamination reagents to monitor the pH throughout the decontamination and buffering procedure, allowing the laboratory technologist to visually see when neutralization has been achieved.

BIBLIOGRAPHY

1. Babakhani, F., Warren, N., Henderson, D., Dalton, H. Effect of Transportation and Acid Neutralization on Recovery of *Mycobacteria* from Processed Specimens. *Am. J. Clin. Pathol.* 1995. 104(1):65-68.
2. Chapmas, J.S., Bernard, J. S. The Tolerance of Unclassified Bacteria. Limits of pH Tolerance. *Am. Rev Respir Dis.* 1962 86:582-583.
3. Kent, P., Kubica, G.P., *Public Health Mycobacteriology: A Guide for the Level III Laboratory.* Centers for Disease Control and Prevention (CDC), 1985.
4. Kubica, G.P., et al. "Sputum Digestion and Decontamination with N-acetyl-L-cysteine-Sodium Hydroxide for Culture of *Mycobacteria*." *Am Rev Respir Dis.* 1963. 87:775-779.
5. Kubica, G.P., et al. "Comments on the Use of the New Mucolytic Agent, N-acetyl-L-cysteine, as a Sputum Digestant for the Isolation of *Mycobacteria*." *Am. Rev. Respir. Dis.* 1964. 89:284-286.
6. Lennett, E.H., et al. *Manual of Clinical Microbiology*, American Society for Microbiology, Third Edition, 1980.
7. Vestal, A.L. *Procedures for the Isolation and Identification of Mycobacteria.* C.D.C., 1975, Atlanta, GA.
8. Yegjian, D., Budd V. "Toxic Effect of Sodium Hydroxide on Tubercle Bacilli." *Am. J. Clinical Pathology*, 1952. 22:456-460.

CONTACT

For Technical Assistance email Technical@AlphaTecSystems.com and for Customer Service, email Sales@AlphaTecSystems.com or call [+1] 800.221.6058 or [+1] 360.260.2779 between 8am and 4pm Monday through Friday, Pacific Time.

WARRANTY

Alpha-Tec Systems, Inc. warrants this product to perform as described in the labeling and literature supplied. Alpha-Tec Systems, Inc. disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall Alpha-Tec Systems, Inc. be liable for any consequential damages arising out of aforesaid express warranty.

TRADEMARKS

NAC-PAC®, NPC-67®, XPR-PLUS®, OxA®, PRB™ and CELL-BOND® are trademarks of Alpha-Tec Systems, Inc., 1311 SE Cardinal Court, Suite 170, Vancouver, WA 98683 USA.

PRODUCT CODES

0004302 NAC-PAC RED (2.0%), 5 x 50 ml, NALC, 5 x 0.25 g	0004307 NAC-PAC RED (2.0%), 5 x 200 ml, NALC, 5 x 1.0 g
0004303 NAC-PAC RED (3.0%), 5 x 50 ml, NALC, 5 x 0.25 g	0004308 NAC-PAC RED (4.0%), 5 x 200 ml, NALC, 5 x 1.0 g
0004304 NAC-PAC RED (4.0%), 5 x 50 ml, NALC, 5 x 0.25 g	0004309 NAC-PAC RED (3.0%), 5 x 200 ml, NALC, 5 x 1.0 g
0004305 NAC-PAC RED (2.5%), 5 x 50 ml, NALC, 5 x 0.25 g	0004311 NAC-PAC RED (3.0%), 200 ml
0004306 NAC-PAC RED (2.5%), 5 x 200 ml, NALC, 5 x 1.0 g	



Manufactured by Alpha-Tec Systems Inc.
1311 SE Cardinal Court, Suite 170, Vancouver, WA 98683 USA.



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



GLOSSARY OF SYMBOLS Update Glossary of Symbols



Lot number / Número de porción / Número de lote / Numerodi lotto / Partienummer / Numéro de sort / Het Aantal van de partij



Catalog number / Número de catálogo / Número de catálogo / Numero di catalogo / Katalogzahl / Numéro de catalogue / Het aantal van de catalogus



For in vitro diagnostic use only / Para el uso diagnóstico in vitro solamente / Para in vitro o uso diagnóstico solamente / Solo per uso diagnostico in vitro / Für nur in vitrodiagnosegebrauch / Pour l'usage diagnostique in vitro seulement / Voor kenmerkend slechts gebruik in vitro



Authorized European representative / Representante europeo autorizado / Representante europeo autorizado / Rappresentante europeo autorizzato / Autorisierter europäischer Repräsentant / Representant européen autorisé / Gemachtigde Europese vertegenwoordiger



Caution, consult accompanying documents / Precaución, consulte los documentos de acompañamiento / Cuidado, consulte originais acompanhando / Attenzione. Consulti i documenti di accompagnamento / Vorsicht, beraten Begleiddokumente / Attention, consultez les documents d'accompagnement / Voorzichtigheid, raadpleeg begeleidende documenten



Store between temperatures indicated / Almacén entre las temperaturas indicadas / Loja entre as temperaturas indicadas / Deposito fra le temperature indicate / Speicher zwischen den Temperaturen angezeigt / Magasin entre les températures indiquées / Opslag tussen vermelde temperature



Consult instructions for use / Consulte las instrucciones para el uso / Consulte instruções para o uso / Consulti le istruzioni per uso / Beraten Sie Anwendungsvorschriften / Consultez les instructions pour l'usage / Raadpleeg instructies voor gebruik



Use by date indicated / Uso por la fecha indicada / Uso pela data indicada / Uso entro la data indicate / Gebrauch bis zum der Dattel angezeigt / Utilisation à la datte indiquée / Gebruik door vermelde datum



Manufacturer / Fabricante / Fabricante / Fornitore / Hersteller / Fabricant / Fabrikant



Contains sufficient for <n> tests / Contiene suficiente para <n> las pruebas / Contem suficiente para <n> testes / Contiene sufficiente per <n> le prove / Enthält genügendes für <n> tests / Contient suffisamment pour <n> des essais / Bevat voldoende voor <n> tests