## TABLE OF CONTENTS

English...	1	Product Codes...	5
French...	3	Glossary of Symbols...	5

## INTENDED USE

ETM® Enteric Transport Media Collection/Transport System is designed for collection, transport, and preservation of patient stool specimens for recovery and examination of enteric bacterial pathogens and toxins (Shiga toxins and *Clostridium difficile* A & B).

## SUMMARY

Bacterial intestinal diseases are diagnosed and treated upon positive recovery and isolation of bacterial pathogens. Different circumstances may delay immediate processing of collected specimens. Refrigeration, incubation, or ambient storage of samples until examination does not provide satisfactory recovery of all enteric pathogens. Therefore, the specimen must be in stasis to guarantee optimal recovery. ETM has been shown to be an excellent transport and preservation medium for *Shigella*, *Salmonella*, *Vibrio*, *Yersinia*, enterotoxigenic *Escherichia coli* (ETEC), *E. coli* O157:H7, and *Campylobacter* as well as toxin assays for *C. difficile* and EIA assays for Shiga toxins [*E. coli* O157:H7] with the EHEC assay. From recent studies, ETM with a pH indicator and 0.16% agar, rather than the traditional 5% agar, maintained *Campylobacter* viability for up to 35 days at 4°C and up to 7 days at 25°C. The efficacy of Buffered Glycerol Saline (BGS) has been demonstrated since its introduction by Teague and Clurman in 1916 and later modified in 1931 by Sachs to support the transport and preservation of *Salmonella* and *Shigella*. However, recent studies have revealed the BGS did not maintain the viability of *Campylobacter jejuni* as long as ETM. For transporting or preserving patient stool specimens until examination for bacteria pathogens, especially *Campylobacter* spp., ETM is recommended. ETM is an isotonic and non-nutritive formulation that prevents overgrowth of commensal flora. The medium is also buffered and a pH indicator is utilized to show any acidic shifts. ETM maintains specimens for enteric pathogen culturing for up to 96 hours at room temperature. ETM may be used for the recovery of *C. difficile* Toxin-A testing when receipt of the specimen into the laboratory is expected to be greater than 2 hours after collection. The specimen should be collected and placed into the ETM medium to yield a 1:3 dilution of specimen-to-medium, then transported to the laboratory at room temperature. Upon receipt, the specimen suspended in ETM should be processed for enteric culture and then held at 2–8°C until tested for *C. difficile* Toxin-A & B.

## FOR IN VITRO DIAGNOSTIC USE ONLY

### STABILITY AND STORAGE

ETM has a shelf life of 24 months from the date of manufacture when stored at room temperature (15–30°C).

### USER QUALITY CONTROL

1. Vials should contain stated volume of clear, red-colored media and be free of precipitates.
2. Any vial with media that is cloudy or yellow should be discarded. The yellow color indicates an acidic pH shift; therefore it should not be accepted for testing. A new specimen should be collected.
3. The total bacteria population of a normal stool specimen must be within +/-20 CFU (Colony Forming Units/ml) for 96 hours.

### SPECIMEN COLLECTION AND PREPARATION

1. Multiple samples can be collected over several days and examined by the laboratory to confirm all findings.
2. Patient specimens should be passed into a dry, clean container such as a bedpan, wax paper, plastic wrap or bag, or Styrofoam plate. The stool must NOT be contaminated with urine or water.
3. Slimy, bloody, or watery stool, if present, should be selected and placed into the vial. Place enough specimen to reach the FILL LINE.

4. Carefully break up the specimen and swirl the contents thoroughly using the spork (attached to the cap). Screw on and tighten the cap without further mixing or shaking.
5. Label each vial with appropriate patient information and deliver the specimen to the laboratory without refrigeration or incubation.

## PROCEDURE

**Materials Provided:** Each vial contains the volume as stated on the label. Patient sets include simple, illustrated, multi-language instructions.

**Materials Not Provided:** Sterile loops or swabs, media.

## SPECIMEN PROCESSING

Due to the hazardous nature of patient specimens, all standard laboratory safety procedures should be observed during all phases of examination and subculturing.

1. Store specimens at room temperature and test culture within 96 hours of collection.
2. Wear disposable gloves at all times.
3. Work in biosafety cabinet at all times.
4. Mix specimen well and remove cap.
5. Place the vial on paper towels pre-moistened with a disinfectant.
6. Inoculate appropriate enrichment and/or selective media using a sterile loop or swab. Good microbiology procedures suggest that bacterial pathogens be processed on both enrichment and selective media.
7. Replace cap and tighten.
8. It is recommended that the vial and all other associated material be autoclaved before disposal.

**NOTE:** If the contents of the specimen vial exhibit a yellow-colored solution, an acid shift has taken place and a new specimen should be collected.

ETM is equivalent to Cary-Blair media and is compatible with the functions, uses, and associated procedures designed for Cary-Blair media.

## BIBLIOGRAPHY

1. Cary, SG, Blair, EB. New transport medium for shipment of clinical specimens, I. fecal specimens. *J Bacteriol.* 1964. 88:96–98.
2. Christian, DL, Ederer, GM. Evaluation of bacteriological transport media. *Am J Med Technol.* 1973. 39(1):12–16.
3. Edwards, PR, Ewing, WH. Identification of Enterobacteriaceae. Burgess Publishing Co., Minneapolis, MN. 3rd ed. 1972. pp. 28, 337–339.
4. Eisenberg, HD, Schoenkecht, FD, Von Graevenitz, A. Culmitech 9-Collection and Processing of Bacterial Specimens. American Society for Microbiology, Washington, D.C. 1979.
5. Ewing, WH. Transport methods for enterobacteriaceae and allied bacteria. *Public Health Lab.* 1971. 29:8–23.
6. Ewing, WH, McWhorter, AC, Montague, TS. Transport media in the detection of *Salmonella typhi* in carriers. *Public Health Lab.* 1966. 24:63–65.
7. Luechtefeld, NW, et al. Evaluation of transportation and storage techniques for isolation of *Campylobacter fetus* subsp. *jejuni* from turkey cecal specimens. *J Clin Micro.* 1981. 13:438–443.
8. Martin, WJ, Washington II, J. "Enterobacteriaceae." *Manual of Clinical Microbiology.* American Society of Microbiology, Washington D.C. 3<sup>rd</sup> ed. 1980. pp. 195–219.
9. Neumann, DA, et al. 1971. Cary-Blair, A transport medium for *Vibrio parahemolyticus*. *Amer J Clin Pathol.* 1971. 57:33–34.
10. Sachs, A. Difficulties associated with the bacteriological diagnosis of bacillary dysentery. *J Royal Army Medical Corps.* 1939. 73:234–239.
11. Teague, O, Clurman, AW. A method of preserving typhoid stools for delayed examination and a cooperative study of the efficiency of each Brilliant-green Agar, Eosin Methylene Blue Agar and Endo Agar for isolation of *Typhoid bacilli* from stools. *J Inf Dis.* 1916.18:635–671.

12. Wang, WL, et al. 1983. Evaluation of transport media for *Campylobacter jejuni* in human fecal specimens. J Clin Micro. 1983.18:803–807.
13. Whittier, S, Clynes, N, Della-Latta, P. Cost effective laboratory practices for the detection of *C. difficile* toxin. Poster Session. 1996. A.S.M. Meeting, New Orleans, LA.

**CONTACT**

Alpha-Tec Systems, Inc. offers a complete parasitology line of collection, transport, and concentration systems as well as QC1™ Quality Control Slides, reagents, and stains. For Technical Assistance email [Technical@AlphaTecSystems.com](mailto:Technical@AlphaTecSystems.com) and for Customer Service, email [Sales@AlphaTecSystems.com](mailto:Sales@AlphaTecSystems.com) or call [+1] 800.221.6058 or [+1] 360.260.2779 between 8 am and 4 pm Monday through Friday, Pacific Time.

**WARRANTY**

Alpha-Tec Systems, Inc. warrants this product to perform as described in the labeling and literature supplied. Alpha-Tec Systems, Inc. disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall Alpha-Tec Systems, Inc. be liable for any consequential damages arising out of aforesaid express warranty.

**TRADEMARKS**

ETM® and QC1™ are trademarks of Alpha-Tec Systems, Inc., 1311 SE Cardinal Court, Suite 170, Vancouver, WA 98683 USA.

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7. Luechtefeld, NW, et al. Evaluation of transportation and storage techniques for isolation of *Campylobacter fetus* subsp. *jejuni* from turkey cecal specimens. J Clin Micro. 1981. 13:438–443.
8. Martin, WJ, Washington II, J. "Enterobacteriaceae." Manual of Clinical Microbiology. American Society of Microbiology, Washington D.C. 3rd ed. 1980. pp. 195–219.
9. Neumann, DA, et al. 1971. Cary-Blair, A transport medium for *Vibrio parahaemolyticus*. Amer J Clin Pathol. 1971. 57:33–34.
10. Sachs, A. Difficulties associated with the bacteriological diagnosis of bacillary dysentery. J Royal Army Medical Corps. 1939. 73:234–239.
11. Teague, O, Clurman, AW. A method of preserving typhoid stools for delayed examination and a cooperative study of the efficiency of each Brilliant-green Agar, Eosin Methylene Blue Agar and Endo Agar for isolation of *Typhoid bacilli* from stools. J Inf Dis. 1916. 18:635–671.
12. Wang, WL, et al. 1983. Evaluation of transport media for *Campylobacter jejuni* in human fecal specimens. J Clin Micro. 1983. 18:803–807.
13. Whittier, S, Clynes, N, Della-Latta, P. Cost effective laboratory practices for the detection of *C. difficile* toxin. Poster Session. 1996. A.S.M. Meeting, New Orleans, LA.

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Manufactured by Alpha-Tec Systems, Inc.  
 1311 SE Cardinal Court, Suite 170  
 Vancouver, WA 98683 USA

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MDSS GmbH