## 5% Oxalic Acid

### TABLE OF CONTENTS

English...	1	Product Codes...	5
French...	3	Glossary of Symbols...	5

### INTENDED USE

5% Oxalic Acid is used for decontamination of clinical respiratory source specimens contaminated with *Pseudomonas* spp.

### SUMMARY

Clinical specimens submitted to the laboratory for the isolation of acid-fast mycobacteria are often contaminated with commensal microbial flora. While most of the contaminating bacteria can be eliminated with a sodium hydroxide and N-acetyl-L-cysteine (NALC) procedure, *Pseudomonas* spp. can survive this decontamination method. In 1930, Corper and Uyei described an oxalic acid decontamination procedure of clinical specimens. In 1993, Whittier et al. demonstrated the effectiveness of the oxalic acid decontamination procedure on *Pseudomonas aeruginosa* in cystic fibrosis patients.

### FOR IN VITRO DIAGNOSTIC USE ONLY

### PRECAUTIONS

- Oxalic Acid is POISONOUS.
- If swallowed, dilute by drinking milk or water. Call a physician or local poison center IMMEDIATELY.
- Avoid contact with skin and eyes. Should contact occur, flush immediately with water. Contact a physician if irritation occurs.
- Observe ALL laboratory standards for the handling of patient specimens.

### STABILITY AND STORAGE

5% Oxalic Acid is stable to the stated expiration date when stored at room temperature (15-30°C).

### USER QUALITY CONTROL

Each bottle of 5% Oxalic Acid should contain the quantity listed on the label and should be clear and free of precipitates.

### SPECIMEN COLLECTION AND PREPARATION

Appropriate specimens for the detection of *Mycobacterium* spp. should be collected according to prescribed standards and delivered to the laboratory in a safe and timely manner. Refer to local procedural guidelines for this information.

### PROCEDURE

**Materials Provided:** 5% Oxalic Acid.

**Materials Not Provided:** Centrifuge, centrifuge tubes (50 ml), microscope slides, sterile pipettes, vortex mixer, NAC-PAC® RED, NPC-67® Neutralizing Buffer or XPR-PLUS® Neutralizing Buffer, PRB™ Pellet Resuspension Buffer, CELL-BOND® Slides, AFB Stains, TB Media.

### SPECIMEN PROCESSING

**\*\*ATTENTION:** Handle all specimens with extreme caution. Possible Pathogens. Wear gloves at ALL TIMES. Work in biosafety cabinet.

- Perform the NALC/NaOH digestion and decontamination procedure on the respiratory specimen. (Call your Alpha-Tec Account Executive for a complete list of our digestion/decontamination and buffering products.)
- If the specimen is identified as having *Pseudomonas* spp. contamination through either a positive contamination control plate or rapid broth test, or patient history of previous *Pseudomonas* infection, proceed with the 5% Oxalic Acid decontamination procedure.

### OXALIC ACID DECONTAMINATION PROCEDURE

- Use up to 3 ml from the original NALC/NaOH processed specimen and add it to a sterile 50ml centrifuge tube.
- Add 3 to 5 ml of 5% Oxalic Acid and vortex for 30 to 60 seconds.
- Allow to stand for 30 minutes, vortexing every 10 minutes.

- Using a sterile pipette, add NAC-PAC RED until the solution turns pink. Re-cap the tube and mix gently. **NOTE:** Any concentration of NAC-PAC RED can be used for this procedure. It is not necessary to add the NALC powder to the NAC-PAC RED.
- Add NPC-67 Neutralizing Buffer or XPR-PLUS Neutralizing Buffer until effective neutralization is indicated by a color change from pink to colorless. Tighten cap, and swirl to mix.
- Centrifuge the specimen tubes at 3000 xg for 15 minutes. It is recommended, but not required, to use a refrigerated centrifuge. (Each laboratory must check the centrifuge head radius and use an appropriate nomogram for proper speed selection (rpm) to achieve the desired relative centrifugal field of 3000 xg.)
- Working in a biosafety hood, pour off all supernatant into a splash-proof container holding an appropriate disinfectant. Use an appropriate disinfectant or flame to disinfect any contamination on the lip of the specimen tube. Do not allow the disinfectant to run down inside the specimen tube.
- Resuspend the pellet with 0.5 to 1.0 ml of PRB Pellet Resuspension Buffer. Do not resuspend the pellet with any other reagents, water, or saline.  
**Option 1:** To maximize time to detection for rapid broth detection systems, resuspend the pellet with 1.0 ml of PRB.  
**Option 2:** Depending on the needs of your laboratory, the pellet may be resuspended with 0.5 ml of PRB to create a more concentrated sample for acid-fast smear staining. Once the smears have been made, add an additional 1 ml of PRB for rapid broth detection methods and media inoculation.
- Mix the sediment and buffer well and inoculate the liquid broth for your rapid broth detection equipment per the manufacturer's instructions.
- Place two drops of the sediment onto the surface of each of the TB media used. A contamination control plate (BAP or TSA) can be inoculated at this point and incubated at 35-37°C for 48 hours.
- Make smears for AFB staining. Use coated CELL-BOND Slides or appropriate sterile albumin adhesive solutions to attach the specimen to the slide. Dry the smears and proceed with AFB staining per the manufacturer's directions. (Call Alpha-Tec Systems, Inc. for a complete list of AFB Stains). **NOTE:** An AFB Quality Control Slide should be stained in conjunction with the patient smears to verify the staining technique and components (QC1™ AFB Slide #0003240).
- Add an additional 1 to 2 ml of PRB to the unused portion of the specimen and refrigerate at 2-8°C for further diagnostic procedures or re-decontamination if the cultures or rapid detection methods indicate that a contaminate is present.

### EXPECTED RESULTS

If Mycobacteria are present in the clinical specimen and processed according to the procedures listed within this document, the recovery of cultivable, viable, and clinically significant mycobacteria can be expected.

### LIMITATIONS OF PROCEDURES

Timing of the decontamination step, proper buffering, speed and timing of the centrifugation step, proper decanting and addition of the Pellet Resuspension Buffer to the pellet are vital to the recovery of mycobacteria. Failure to follow the listed procedures may result in decreased numbers of mycobacteria or total loss of mycobacteria resulting in an inaccurate culture report.

### BIBLIOGRAPHY

- CDC. Verbal Consultation, May 18, 1989.
- Corper, H.J., and N. Uyei. "Oxalic acid as a reagent for isolating tubercule bacilli and a study of the growth of acid-fast nonpathogens on different mediums with their reactions to chemical reagents." J Lab Clin Med. 1930.15:348-369.
- Lenette, E.H. et al. Manual of Clinical Microbiology. American Society of Microbiology, Washington, D.C. 1980. pp. 150-179.

## 5% Oxalic Acid

4. Vestal, A. Procedures for the isolation and identification of mycobacteria. U.S. Dept. of Health. CDC Publication No. 79-8230, Center for Disease Control, Atlanta, GA. 1975. pp.21-31.
5. Whittier, S., et al. "Improved recovery of mycobacteria from respiratory secretions of patients with cystic fibrosis". J Clin Microbiol. 1993. 31:861-864.

### CONTACT

Alpha-Tec Systems, Inc. offers a complete line of reagents, stains, and QC1™ Quality Control Slides for AFB, Parasitology, Bacteriology, and Mycology processing, as well as O&P collection systems and concentration devices for Parasitology. For Technical Assistance, email [Technical@AlphaTecSystems.com](mailto:Technical@AlphaTecSystems.com), and for Customer Service, email [Sales@AlphaTecSystems.com](mailto:Sales@AlphaTecSystems.com), or call either [+1] 800.221.6058 (USA) or [+1] 360.260.2779 between 8AM and 4PM Monday through Friday, Pacific Time.

### WARRANTY

This product is warranted by Alpha-Tec Systems, Inc. to perform as described in the labeling and literature supplied. Alpha-Tec Systems, Inc. disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall Alpha-Tec Systems, Inc. be liable for any consequential damages arising out of aforesaid express warranty.

### TRADEMARKS

CELL-BOND®, NAC-PAC®, NPC-67®, PRB™, QC1™ and XPR-PLUS® are trademarks of Alpha-Tec Systems, Inc., 1311 SE Cardinal Court, Suite 170, Vancouver, WA 98683 USA.

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**BIBLIOGRAPHIE**

1. CDC, Verbal Consultation, May 18, 1989.
2. Corper, H.J., and N. Uyei. "Oxalic acid as a reagent for isolating tubercule bacilli and a study of the growth of acid-fast nonpathogens on different mediums with their reactions to chemical reagents." J Lab Clin Med. 1930.15:348-369.
3. Lennette, E.H. et al. Manual of Clinical Microbiology. American Society of Microbiology, Washington, D.C. 1980. pp. 150-179.
4. Vestal, A. Procedures for the isolation and identification of mycobacteria. U.S. Dept. of Health. CDC Publication No. 79-8230, Center for Disease Control, Atlanta, GA. 1975. pp.21-31.
5. Whittier, S., et al. "Improved recovery of mycobacteria from respiratory secretions of patients with cystic fibrosis". J Clin Microbiol. 1993. 31:861-864.

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Manufactured by Alpha-Tec Systems, Inc.  
1311 SE Cardinal Court, Suite 170  
Vancouver, WA 98683 USA

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MDSS GmbH