TABLE OF CONTENTS

English...	1	Product Codes...	5
French...	3	Glossary of Symbols...	5

INTENDED USE

MUCOSOL™ is used in qualitative procedures as a mucolytic agent in the digestion of disulfide bonds in mucous for qualitative aerobic, mycological and mycobacterial culture, identification procedures, and nucleic acid extraction methods.

SUMMARY

The Dithiothreitol technique was described in 1963 by W.W. Cleland. Reep and Kaplan evaluated Dithiothreitol and N-acetyl-L-cysteine as mucolytic agents for sputum specimens and demonstrated that either agent was successful in liquefying and decontaminating sputa for acid-fast culture. They reported that by using sodium citrate without sodium hydroxide as the diluent, the toxicity to fungus was eliminated. Dithiothreitol is also an effective agent in liquefying sputa for molecular diagnostics.

FOR IN VITRO DIAGNOSTIC USE ONLY
PRECAUTIONS

This product should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers and media after use. Directions should be read and followed carefully.

STABILITY AND STORAGE

MUCOSOL can be stored in its original container at room temperature. When diluted, refrigerate at 2-8°C. The diluted solution should be used immediately. Bring product to room temperature prior to use. If proper sterile technique is used, the diluted solution is stable for at least 48 hours at 2-8°C. If product freezes, a crystal precipitate may form, but will not affect product performance. Place in a water bath to dissolve crystals.

QUALITY CONTROL

Any product showing cloudiness, turbidity, precipitation, or discoloration should be discarded. Quality controlled microorganisms should be utilized to verify procedures, media and reagents as appropriate for your regulatory agency or local procedural guidelines.

SPECIMEN COLLECTION AND PREPARATION

Specimens should be collected and transported to the laboratory without delay and protected from excessive heat or cold. Specimens for AFB should be collected without any fixatives.

PROCEDURE

Materials Provided: MUCOSOL (Dithiothreitol), 10 ml

Materials Not Provided: Centrifuge, centrifuge tubes, Vortex mixer, pipettes, microscope slides, culture media, nucleic acid extraction reagents, AFB processing reagents, Bovine Albumin.

For Routine Cultures and Nucleic Acid Extraction

1. Add the contents of one MUCOSOL vial to 90 ml of sterile distilled water and mix gently. Or, dilute 1:10 to an appropriate volume for the application.
2. Line up specimens and centrifuge tubes in a laminar flow hood.
3. Loosen specimen container caps. Work in sets equivalent to a centrifuge load.
4. Transfer the specimen to a sterile 50 ml centrifuge tube or other appropriate volume of centrifuge tube. Do not use more than 10 ml of specimen. NOTE: The amount of specimen to be processed in the centrifuge tube must never exceed more than 1/5 the volume

of the tube (i.e., for a 50 ml centrifuge tube, 10 ml of specimen is the maximum.). If it is necessary to process the entire specimen, use more than one centrifuge tube and combine the sediments after centrifugation.

5. Add a volume of the MUCOSOL to the specimen equal to the size of the specimen (i.e., if the specimen is approximately 5 ml, add 5 ml of the MUCOSOL solution to the same specimen tube).
6. Tighten the caps on the centrifuge tubes. Mix each specimen on a Vortex until liquefied (5 to 20 seconds per specimen).
7. Allow each specimen to stand for 15-30 minutes.
8. Centrifuge at 3000 xg for 15 minutes. (Each laboratory must check the centrifuge head radius, and then use an appropriate nomogram for proper speed selection (rpm) to achieve the desired relative centrifugal field of 3000 xg.) It is recommended but not required to use a refrigerated centrifuge.
9. Working in a biosafety hood, pour off half of the supernatant. The resulting volume is now equal to the starting volume of the sputum. Mix with a vortex to resuspend the pellet.
10. Process the sample as liquefied sputum for routine culture or molecular extraction by the methods approved by your laboratory.

For Acid-Fast Cultures Only

1. Add the contents of one MUCOSOL vial to 90 ml of sterile distilled water and mix gently.
2. Line up specimens and centrifuge tubes in a laminar flow hood.
3. Loosen specimen container caps. Work in sets equivalent to a centrifuge load.
4. Transfer the specimen to a sterile 50 ml centrifuge tube. Do not use more than 10 ml of specimen. NOTE: The amount of specimen to be processed in the centrifuge tube must never exceed more than 1/5 the volume of the tube (i.e., for a 50 ml centrifuge tube, 10 ml of specimen is the maximum). If it is necessary to process the entire specimen, use more than one centrifuge tube and combine the sediments after centrifugation.
5. Add a volume of the MUCOSOL to the specimen equal to the size of the specimen (i.e., if the specimen is approximately 5 ml, add 5 ml of the MUCOSOL solution to the same specimen tube).
6. Tighten the caps on the centrifuge tubes. Mix each specimen on a Vortex until liquefied (5 to 20 seconds per specimen).
7. Allow each specimen to stand for 15 minutes.
8. Centrifuge at 3000 xg for 15 minutes. (Each laboratory must check the centrifuge head radius, and then use an appropriate nomogram for proper speed selection (rpm) to achieve the desired relative centrifugal field of 3000 xg.) It is recommended but not required to use a refrigerated centrifuge.
9. Working in a biosafety hood, pour off all the supernatant into a splash-proof container. The container should contain an appropriate disinfectant.
10. To the sediment, add an equal amount of an appropriate sodium hydroxide digestion solution.
11. Tighten the caps on the centrifuge tubes. Mix each specimen on a vortex until liquefied (5 to 20 seconds per specimen).
12. Allow each specimen to stand for 15 minutes.
13. Neutralize the sodium hydroxide with an appropriate buffer solution. For best results, ensure the pH of the solution is neutral. Tighten the cap and swirl by hand to mix.
14. Centrifuge at 3000 xg for 15 minutes.
15. Within a hood, pour off all the supernatant into a splash-proof container. The container should contain an appropriate disinfectant. Using a sterile pipette, add 1 or 2 ml of the 0.2% Bovine Albumin Fraction V Solution and mix to resuspend the sediment. Or, resuspend the sediment in an appropriate buffer and use a coated slide.
16. Inoculate the sediment onto the surface of each of the TB media used.
17. Make a smear for acid-fast staining. Dry the smears and proceed with appropriate acid-fast staining per the manufacturer's instructions.

For Mycological Cultures

1. Add the contents of one MUCOSOL vial to 90 ml of 1.47% Sodium Citrate (Trisodium Citrate) Solution.
2. Line up specimens and centrifuge tubes in laminar flow hood.
3. Loosen specimen container caps. Work in sets equivalent to a centrifuge load.
4. Transfer the specimen to a sterile 50 ml centrifuge tube. Do not use more than 10 ml of specimen. NOTE: The amount of the specimen to be processed in the centrifuge tube must never exceed more than 1/5 the volume of the tube (i.e., for a 50 ml centrifuge tube, 10 ml of specimen is the maximum). If it is necessary to process the entire specimen, use more than one centrifuge tube and combine the sediments after centrifugation.
5. Add the volume of the MUCOSOL to the specimen equal to the size of the specimen (i.e., if the specimen is approximately 5 ml, add 5 ml of the MUCOSOL Solution to the same specimen tube).
6. Tighten the caps on the centrifuge tubes. Mix each specimen on a Vortex until liquefied (5 to 20 seconds per specimen).
7. Allow each specimen to stand for 15 minutes.
8. Centrifuge at 3000 xg for 15 minutes. (Each laboratory must check the centrifuge head radius, and then use an appropriate nomogram for proper speed selection (rpm) to achieve the desired relative centrifugal field of 3000 xg.) It is recommended but not required to use a refrigerated centrifuge.
9. Working in a biosafety hood, pour off all but approximately 1-2 ml of the supernatant into a splash-proof container. The container should contain an appropriate disinfectant.
10. Using a sterile pipette, resuspend the sediment.
11. Make all appropriate smears and inoculate an appropriate quantity of the sediment to mycological media.

EXPECTED RESULTS

If aerobic bacteria, mycobacteria, fungi, or viral material are present in the clinical specimen and processed according to the procedures listed with in this document, the recovery of cultivatable, viable, and clinically significant organisms and nucleic acids can be expected.

LIMITATIONS OF PROCEDURES

This product is only part of the overall scheme for identification. Procedures for biochemical and serological tests for identification may be found in appropriate references.

BIBLIOGRAPHY

1. Cleland W.W., "Dithiothreitol: A New Protective Reagent for SH Groups", BIOCHEMISTRY, Vol 3 No. 4, April 1964.
2. Hirsch, S. Roger, Joyce E. Zastrow and Ross C. Kory, "Sputum Liquefying Agents: A Comparative In Vitro Evaluation", J. LAB. CLIN. MICRO., Aug 1969.
3. Kubica, G.P., et al. Sputum Digestion and Decontamination with N-acetyl-L-cysteine-Sodium Hydroxide for Culture of Mycobacteria. Am Rev Respir Dis., 87:775-779, May 1963.
4. Kubica, G.P. et al. Comments on the Use of the New Mucolytic Agent, N-acetyl-L-cysteine, as a Sputum Digestant for the Isolation of Mycobacteria. Am Rev Respir Dis., 89:284-286, February 1964.
5. Lennette, E.H., A. Balows, W.J. Hausler, Jr. and J.P. Truant, MANUAL OF CLINICAL MICROBIOLOGY, 3rd edition, 1980.
6. Reep, B.R. and W. Kaplan, "The Effect of Newer Tubercle Bacillus Digestion and Decontamination Procedures on Fungi Causing Pulmonary Disease", MYCOPATHOLOGIA AND MYCOLOGIA APPLICATA, Vol. 46: 325-334, 1972.
7. Reep, B.R. and W. Kaplan, "The Use of N-acetyl-L-cysteine and Dithiothreitol to Process Sputa for Mycological and Fluorescent Antibody Examination", H.L.S., Vol. 9, No. 2, 1972.
8. Vestal, Annie L., "Procedures for the Isolation and Identification of Mycobacteria", C.D.C., Atlanta, 1975.
9. "Processing of Sputum Specimens for Nucleic Acid Extraction", CDC, 2020.

CONTACT

CalibreScientific US, Inc. offers a complete line of reagents, stains, and QC1™ Quality Control Slides for AFB, Parasitology, Bacteriology, and Mycology processing, as well as O&P collection systems and concentration devices for Parasitology. For Technical Assistance, email Technical@AlphaTecSystems.com, and for Customer Service, email Sales@AlphaTecSystems.com, or call either [+1] 800.221.6058 (USA) or [+1] 360.260.2779 between 8AM and 4PM Monday through Friday, Pacific Time.

WARRANTY

This product is warranted by CalibreScientific US, Inc. to perform as described in the labeling and literature supplied. CalibreScientific US, Inc. disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall CalibreScientific US, Inc. be liable for any consequential damages arising out of aforesaid express warranty.

TRADEMARKS:

MUCOSOL™ and QC1™ are trademarks of CalibreScientific AMER IP LLC., 6201 Trust Dr, Holland, OH 43528.

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13. Kubica, G.P. et al. Comments on the Use of the New Mucolytic Agent, N-acetyl-L-cysteine, as a Sputum Digestant for the Isolation of Mycobacteria. *Am Rev Respir Dis.*, 89:284-286, February 1964.
14. Lennette, E.H., A. Balows, W.J. Hausler, Jr. and J.P. Truant, *MANUAL OF CLINICAL MICROBIOLOGY*, 3rd edition, 1980.
15. Reep, B.R. and W. Kaplan, "The Effect of Newer Tubercle Bacillus Digestion and Decontamination Procedures on Fungi Causing Pulmonary Disease", *MYCOPATHOLOGIA AND MYCOLOGIA APPLICATA*, Vol. 46: 325-334, 1972.
16. Reep, B.R. and W. Kaplan, "The Use of N-acetyl-L-cysteine and Dithiothreitol to Process Sputa for Mycological and Fluorescent Antibody Examination", *H.L.S.*, Vol. 9, No. 2, 1972.
17. Vestal, Annie L., "Procedures for the Isolation and Identification of Mycobacteria", C.D.C., Atlanta, 1975.
18. "Processing of Sputum Specimens for Nucleic Acid Extraction", CDC, 2020.

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PRODUCT CODES:

0003510 MUCOSOL (Cleland's Reagent), 15 x 10 ml

 Manufactured by CalibreScientific US, Inc.
 1311 SE Cardinal Court, Suite 170
 Vancouver, WA 98683 USA

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MUCOSOL™ et QC1™ sont des marques déposées de CalibreScientific AMER IP LLC., 6201
MDSS GmbH