**Kinyoun Stain Set with Methylene Blue,  
Kinyoun Carbolfochsin Stain, TB Decolorizer,  
Methylene Blue Counterstain**

**TABLE OF CONTENTS**

English...	1	German...	9
French...	3	Product Codes...	11
Spanish...	5	Glossary of Symbols...	11
Italian...	7		

**INTENDED USE**

Alpha-Tec Systems, Inc. Kinyoun Stain Set with Methylene Blue and components are recommended for use in qualitative procedures for the microscopic detection of mycobacteria.

**SUMMARY**

One of the earliest methods devised for the detection of the tubercle bacillus is the microscopic staining technique. This method is valuable for the rapid initial step in the presumptive diagnosis of mycobacterial infections and provides information about the number of acid-fast bacilli present. Species of *Nocardia* and *Rhodococcus* are also partially acid-fast. Three types of staining procedures are used for the detection of acid-fast organisms: fluorochrome, Ziehl-Neelsen (hot stain), and Kinyoun (cold stain).

The basic Kinyoun Carbolfochsin dye used in this stain interacts with the mycolic acids and waxes in the mycobacterial cell wall, which then fails to wash out with mild acid decolorization. Typical acid-fast bacteria, which retain the dye, appear pink to red. The Methylene Blue Counterstain imparts a contrasting color to detect non acid-fast bacteria in the smear.

**FOR IN VITRO DIAGNOSTIC USE ONLY**
**PRECAUTIONS**

Kinyoun Carbolfochsin Stain is poisonous and may be harmful or fatal if swallowed. May cause irritation to skin, eyes and respiratory tract. Combustible; keep away from heat and flame. Suspected cancer hazard; contains material which may cause cancer. Avoid breathing vapor and eye/skin contact. TB Decolorizer may be harmful or fatal if swallowed. May cause irritation to skin, eyes and respiratory tract. Flammable, keep away from heat and flame. Avoid breathing vapor and eye/skin contact.

**STORAGE**

This product is ready to use, and no further preparation is necessary. Store in original containers at 15-30°C until use.

**USER QUALITY CONTROL**

This product should not be used if the color has changed, the expiration date has passed, or there are other signs of deterioration.

**SPECIMEN COLLECTION AND PREPARATION**

Appropriate specimens for the detection of *Mycobacterium* spp. should be collected according to prescribed standards and delivered to the laboratory in a safe and timely manner. Refer to local procedural guidelines for this information.

**PROCEDURE**

**Materials Provided:** Kinyoun Carbolfochsin Stain, TB Decolorizer, Methylene Blue Counterstain. Each component is also available individually.

**Materials Not Provided:** Pipettes, glass slides, Bunsen burner or slide warmer, microscope, immersion oil, deionized or distilled water, QC1™ AFB Slide (#0003240).

**SPECIMEN PROCESSING**

1. Make a thin smear of the specimen over an area of 2-3 sq cm. Heat fix by warming on a slide warmer (65-75°C) for 10 minutes, or until the smear is dry. If a Bunsen burner is used, pass the slide slowly 3 times over its cone of heat. Do not scorch.
2. Flood the smear with Kinyoun Carbolfochsin Stain for 5 minutes.

3. Rinse with deionized or distilled water and drain.
4. Decolorize with TB Decolorizer for 3 minutes.
5. Rinse with deionized or distilled water and drain.
6. Repeat decolorization for 1-2 minutes or until no red stain appears.
7. Rinse with deionized or distilled water and drain.
8. Flood smear with Methylene Blue Counterstain for 3-4 minutes.
9. Rinse with deionized distilled water and allow to air-dry.
10. Examine microscopically using low power (10x-40x) and confirm at 100x under oil immersion for the presence of acid-fast bacilli.

**EXPECTED RESULTS**

**Positive Reaction:** Acid-fast positive organisms stain pink to red. Mycobacteria appear red, small, slightly curved, possibly beaded or banded, with tapered ends.

**Negative Reaction:** Non acid-fast organisms stain blue.

**LIMITATIONS OF PROCEDURES**

1. Overheating (burning) during fixation should be avoided.
2. Do not stain smears that have only been air-dried. Smears must also be "fixed".
3. Smears should not be too thick. After air-drying, examine under a microscope.
4. After staining, it is essential that the back surface of the slide is wiped clean.
5. If washing with distilled water is not done adequately, crystallization of the stain may appear on the slide.
6. Scan at least 100 fields before reporting as negative.
7. Some non-tuberculous species, such as *Mycobacterium avium* complex, appear pleomorphic, usually coccoid. Positive findings should be based only on typical forms but atypical cells should be noted.
8. Atypical rods may represent other partially acid-fast organisms such as species of *Nocardia* or *Rhodococcus*. A weaker acid or a shorter decolorizing period is used to detect these organisms.
9. The sensitivity of the direct acid-fast smear examination for the diagnosis of mycobacterial infection is lower than that of culture methods. Therefore, cultures should be performed on all specimens.

**SPECIFIC PERFORMANCE CHARACTERISTICS**

This product has been tested using the QC1™ AFB Quality Control Slide and has been found to yield acceptable stain results. Appropriate testing should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

**BIBLIOGRAPHY**

1. Baron, E.J., L.R. Peterson, and S.M. Finegold. 1994. Bailey and Scott's Diagnostic Microbiology. 9<sup>th</sup> Ed. Mosby, St. Louis, MO.
2. Bennedson, J. and S.O. Larsen. 1966. Scand. J. Respir. Dis. 47:114-120.
3. Bishop, P.J. and G. Neumann. 1970. TUBERCLE. 51:196-206.
4. Garcia. S. and D.A. Bruckner. 1997. Diagnostic Medical Parasitology. 3<sup>rd</sup> Ed. ASM, Washington, D.C.
5. Kinyoun, J.J. 1915. Am. J. Public Health. 5:867-870.
6. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover. 1995. Manual of Clinical Microbiology. 6<sup>th</sup> Ed. ASM, Washington, D.C.

**CONTACT**

Alpha-Tec Systems, Inc. offers a complete line of reagents, stains, and QC1™ Quality Control Slides for AFB, Parasitology, Bacteriology, and Mycology processing, as well as O&P collection systems and concentration devices for Parasitology. For Technical Assistance, email [Technical@AlphaTecSystems.com](mailto:Technical@AlphaTecSystems.com), and for Customer Service, email [Sales@AlphaTecSystems.com](mailto:Sales@AlphaTecSystems.com), or call either [+1] 800.221.6058 (USA) or [+1] 360.260.2779 between 8AM and 4PM Monday through Friday, Pacific Time.

[+1] 800.221.6058 United States  
[+1] 360.260.2779 International

Directions For Use for the following:

**Kinyoun Stain Set with Methylene Blue,  
Kinyoun Carbolfuchsin Stain, TB Decolorizer,  
Methylene Blue Counterstain**

**WARRANTY**

This product is warranted by Alpha-Tec Systems, Inc. to perform as described in the labeling and literature supplied. Alpha-Tec Systems, Inc. disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall Alpha-Tec Systems, Inc. be liable for any consequential damages arising out of aforesaid express warranty.

**TRADEMARKS**

QC1™ is a trademark of Alpha-Tec Systems, Inc., 1311 SE Cardinal Court, Suite 170, Vancouver, WA 98683 USA.

**Kinyoun Stain Set with Methylene Blue,  
Kinyoun Carbolfuchsin Stain, TB Decolorizer,  
Methylene Blue Counterstain**

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**Kinyoun Stain Set with Methylene Blue,  
Kinyoun Carbol-fuchsin Stain, TB Decolorizer,  
Methylene Blue Counterstain****LITERATUR**

1. Baron, E.J., L.R. Peterson, and S.M. Finegold. 1994. Bailey and Scott's Diagnostic Microbiology. 9th Ed. Mosby, St. Louis, MO.
2. Bennedson, J. and S.O. Larsen. 1966. Scand. J. Respir. Dis. 47:114-120.
3. Bishop, P.J. and G. Neumann. 1970. TUBERCLE. 51:196-206.
4. Garcia. S. and D.A. Bruckner. 1997. Diagnostic Medical Parasitology. 3rd Ed. ASM, Washington, D.C.
5. Kinyoun, J.J. 1915. Am. J. Public Health. 5:867-870.
6. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover. 1995. Manual of Clinical Microbiology. 6TH Ed. ASM, Washington, D.C.

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**PRODUCT CODES**

0003524	Kinyoun Carbolfuchsin Stain, 250 ml
0003524C	Kinyoun Carbolfuchsin Stain, 8 x 250 ml
0003525	TB Decolorizer, 250 ml
0003525C	TB Decolorizer, 8 x 250 ml
0003526	Methylene Blue Counterstain, 250 ml
0003526C	Methylene Blue Counterstain, 8 x 250 ml
0003566	Kinyoun Stain Set with Methylene Blue, 1 set, 3 x 250 ml

Manufactured by Alpha-Tec Systems, Inc.  
1311 SE Cardinal Court, Suite 170  
Vancouver, WA 98683 USA

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QC1™ ist ein Warenzeichen von Alpha-Tec Systems, Inc., 1311 SE Cardinal Court, Suite 170, Vancouver, WA 98683 USA.
MDSS GmbH