TABLE OF CONTENTS

English...	1	Product Codes...	6
Italian...	4	Glossary of Symbols...	6

INTENDED USE

Parasite Recovery System (PRS™) is a complete system for concentration of helminth eggs, larvae and protozoan cysts. Ethyl Acetate and Triton X-100™ are reagents used in the concentration procedure.

SUMMARY

The diagnosis of intestinal parasitic infections requires appropriate collection, transport, concentration, staining and identification of parasites from fecal specimens. A wide variety of sedimentation, concentration, and flotation methods have been described. The formalin-ether sedimentation technique (Ritchie, et al.) has been modified to avoid the flammability, storage and disposal of ether. The modified procedure replaces the ether with ethyl acetate. This procedure can be performed on specimens that have been fixed in Buffered Formalin, SAF, PROTO-FIX®, or on fresh specimens where the recommended fixative has been added during the processing.

The Parasite Recovery System utilizes existing methodologies and standardized components to accommodate specimen preparation for subsequent microscopic examination.

1. The addition of the surfactant to preserved specimens can help reduce the adhesive forces of the mucus, break up the lumps within the sample and increase settling of the parasitic eggs.
2. The PRS funnel/filtration unit has a precision molded filtration screen, which allows helminth eggs, larvae, and protozoan trophozoites and cysts to pass through the screen but retains the larger diameter particulate matter on the top. Most of the macroscopic fecal debris will be discarded with the funnel/filtration unit.
3. Ethyl Acetate is less flammable than ethyl ether and as such, is recommended as a replacement to ether. The Ethyl Acetate will dissolve the fat and aid in the separation of the fecal debris from a concentrated sediment containing parasites, if present in the sample.

FOR IN VITRO DIAGNOSTIC USE ONLY

PRECAUTIONS

1. **CAUTION!** Ethyl Acetate is **FLAMMABLE**. Perform all procedures in a well-ventilated area. Do not allow any open flames or ignition devices in the room when these procedures are being performed. Avoid prolonged breathing of fumes. Avoid skin contact or contact with eyes. Wear gloves and eye protection at all times while performing these procedures.
2. Observe standard Good Laboratory Practices in handling and disposing of bio-hazardous clinical specimens and laboratory reagents. Refer to your facility safety director for specific details.
3. Do not use the product if the expiration date on the reagents has been exceeded.

SPECIMEN COLLECTION & PREPARATION

1. Specimens preserved in Buffered Formalin, SAF, PROTO-FIX, or fresh samples may be processed with the PRS. The specimen must be fixed for a minimum of 60 minutes to assure adequate fixation of the sample. The specimen should be stored at room temperature.
2. An appropriate clinical patient sample, collected, preserved/fixed, and transported properly is important to the recovery of helminth eggs, larvae and protozoan cysts. Please refer to appropriate references for the collection and transport methods.
3. Always mix the sample well.
4. The appropriate volume of sample is 2-3 grams of fecal matter in 13-15 ml of fixative.

PROCEDURE

Materials Provided: PRS Funnels, and/or Triton X-100, and/or ethyl acetate, and/or centrifuge tubes and caps.

Materials Not Provided: Physiological saline, microscope, microscope slides, coverslips, centrifuge, cotton-tipped applicator sticks, pipettes, 10% Buffered Formalin.

SPECIMEN PROCESSING

PROTO-FIX Preserved Specimens

Refer to the PROTO-FIX Directions For Use for recommended processing procedures.

Specimen Processing CONSED® Concentrating Solution

Refer to the CONSED Concentrating Solution Directions For Use for recommended processing procedures.

Ethyl Acetate Method

Unpreserved Specimens

For optimum results, it is recommended that specimens be preserved at the time of collection. Unpreserved specimens delayed in transport may have limited diagnostic value.

1. Transfer 3-5 grams of unpreserved stool into 15 ml preservative of choice. Mix stool with preservative thoroughly and break up any lumps or fecal masses. The stool/preservative mixture should stand for a minimum of 60 minutes for adequate fixation. Follow procedure for use below for the specific fixative.

Immediate Processing of Unpreserved Specimens

1. Transfer 5-6 grams of unpreserved stool into 10-15 ml of physiological saline. Mix stool/saline mixture thoroughly, breaking up any lumps or fecal masses.
2. Add 4 drops of Triton X-100 to the mixture. (Up to 8 drops of Triton X-100 may be added if the specimen is highly mucoid.)
3. Cap and mix the contents thoroughly by shaking.
4. Insert one of the PRS filtering devices into the top of the 15 ml centrifuge tube.
5. Pour fecal suspension through the filtering device into the centrifuge tube. Usually 3 ml in the centrifuge tube is sufficient unless the fecal suspension is thin.
6. Discard filtering device, add 10 ml physiological saline and centrifuge for 10 minutes at 500-600 xg (not rpm). Decant the supernatant fluids, retaining the sediment. About 1 ml of sediment should be present.
7. Resuspend the sediment in 10 ml buffered formalin. Allow mixture to stand for at least 5 minutes before proceeding. (At this point, the mixture in the centrifuge tube may be capped and saved until a later time.)
8. Add approximately 3 ml of Ethyl Acetate. Cap the tube and shake for 30 seconds. Invert the tube while shaking. **CAUTION:** Pressure may build up within the tube during shaking. Carefully release the pressure by opening the cap on the centrifuge tube away from you.
9. Centrifuge the tube for 10 minutes at 500-600 xg. Examination of the tube after centrifugation should reveal 4 distinct layers from the top down: (1) a layer consisting of Ethyl Acetate, (2) a "plug" of fecal debris, (3) a discolored aqueous layer, (4) a sediment layer, containing the parasites.
10. The final sediment remaining should be approximately 0.25 ml. Hold the tube in a vertical position. Free the plug of debris by ringing with a wooden applicator stick. Decant the upper layers, leaving the sediment. Do not turn the tube upright until the sides of the tube have been leaned with cotton-tipped swabs.
11. Transfer a portion of the sediment to a clean glass microscope slide and prepare the mount of choice. Examine microscopically. A portion of the sediment may be used for detection of *Cryptosporidium*. Consult an appropriate reference for proper preparation and examination.

Formalin Preserved Specimens

1. Thoroughly mix the formalized specimen.
2. Add 5 drops of the Triton X-100 reagent to the specimen in the formalin vial.
3. Recap the fixative vial. Mix the fixed specimen with the surfactant thoroughly by shaking 10 to 20 seconds.
4. Appropriately label one side of the 15 ml centrifuge tube. Place the labeled centrifuge tube in a rack and insert one of the PRS funnel/filtration units into the tip of the 15 ml centrifuge tube.
5. Carefully pour the fecal suspension through the filtration device into the centrifuge tube. A minimum of 1 ml of filtered sediment is necessary to achieve a quality specimen. Approximately 3 ml of sediment is usually adequate and up to 10 ml of filtered sediment may be necessary for watery fecal samples. Do not force the fecal/fixative solution through the PRS filtration device.
6. Discard the PRS filtration unit into appropriate disposal container.
7. Add 10 ml of physiological saline to the filtered sample (in the 15 ml centrifuge tube) and centrifuge 500 – 600 xg for 10 minutes.
8. Decant the supernatant, retaining the sediment.
9. Resuspend the sediment to the 9 ml mark on the centrifuge tube with 10% Buffered Formalin and mix the specimen well.
10. Add 4 ml of Ethyl Acetate or substitute. Cap the tube, invert, and shake vigorously for 30 seconds. Remove the cap with care, pointing the tube away from you, since pressure can build up within the tube during shaking.
11. Centrifuge at 500-600 xg for 10 minutes. Following centrifugation, four (4) layers should result (from the top down): (1) Ethyl Acetate, (2) a "plug" of fecal debris and fat, (3) Formalin, (4) Sediment layer.
12. Holding the centrifuge tube in a vertical position, free the plug of debris from the sides of the tube by ringing the tube with a wood applicator stick. Carefully pour the top three layers into an appropriate discard container. While the tube is still tipped in the decanting position, use a cotton-tipped swab to remove debris from the sides of the tube. Do not turn the tube upright until the sides of the tube have been thoroughly cleaned. This step is very important. If the sides are not cleaned thoroughly with the swab, lipid droplets can mix with the sediment, making microscopic examination more difficult.
13. Prepare wet mounts of the sediment on a clean glass microscope slide. Examination should be performed within 30 minutes. If mounts are to be prepared later, a small amount of 10% Buffered Formalin may be added to the sediment tube and tightly capped. Consult appropriate references for the proper examination of the sediment and identification of parasites.

SAF Preserved Specimens

1. Thoroughly mix the preserved specimen.
2. Add 5 drops of the Triton X-100 reagent to the specimen in the SAF vial.
3. Recap the fixative vial. Mix the fixed specimen with the surfactant thoroughly by shaking 10-20 seconds.
4. Appropriately label one of the 15 ml centrifuge tubes. Place the labeled centrifuge tube in a rack and insert one of the PRS funnel/filtration units into the top of the 15 ml centrifuge tube.
5. Carefully pour the fecal suspension through the filtration device into the centrifuge tube. Approximately 3 ml of sediment is usually adequate and up to 10 ml of filtered sediment is necessary for watery fecal samples. Do not force the fecal/fixative solution through the PRS filtration device.
6. Discard the PRS filtration unit into an appropriate disposal container.
7. Add 10 ml of physiological saline to the filtered sample (in the 15 ml centrifuge tube) and centrifuge between 500 - 600 xg for 10 minutes.
8. Decant the supernatant, retaining the sediment.
9. If permanent stained slides are necessary, prepare as follows: (1) Transfer some of the sediment onto a CELL-BOND® Slide or onto a clean glass microscope slide and mix with a drop of Mayer's Albumin (#0003355). (2) Mix the specimen well and smear a film

over the slide. Consult appropriate references for proper staining techniques.

10. Resuspend the sediment to the 9 ml mark on the centrifuge tube with 10% Buffered Formalin and mix the specimen well.
11. Add 4 ml of Ethyl Acetate. Cap the tube, invert the tube and shake vigorously for 30 seconds. Remove the cap with care, pointing the tube away from you, since pressure can build up within the tube during shaking.
12. Centrifuge at 500-600xg for 10 minutes. Following Centrifugation four (4) layers should result (from the top down): (1) Ethyl Acetate, (2) A "plug" of fecal debris, and fat, (3) Formalin, (4) Sediment layer.
13. Holding the centrifuge tube in a vertical position, free the plug of debris from the sides of the tube by ringing the tube with a wood applicator stick. Carefully pour the top three layers into an appropriate discard container. While the tube is still tipped in the decanting position, use a cotton-tipped swab to remove debris from the sides of the tube. Do not turn the tube upright until the sides of the tube have been thoroughly cleaned. This step is very important. If the sides are not cleaned thoroughly with the swab, lipid droplets can mix with the sediment, making microscopic examination more difficult.
14. Prepare wet mounts of the sediment on a clean glass microscope slide. Examination should be performed within 30 minutes. If mounts are to be prepared later, a small amount of 10% Buffered Formalin may be added to the sediment tube and tightly capped. Consult appropriate references for the proper examination of the sediment and identification of parasites.

SPECIFIC PERFORMANCE CHARACTERISTICS

The Parasite Recovery System (PRS) should provide comparable results to the formalin-ethyl acetate concentration procedure described by Ritchie and as modified by Young, et al.

Notes:

1. If the fresh fecal sample is watery, the specimen should first be centrifuged for 3-5 minutes at 450-500 xg. Carefully decant the supernatant and use the sediment for concentration procedures.
2. The PRS Funnel/Filtration Unit has a precision molded filtration screen, which allows helminth eggs and larvae, and protozoan trophozoites and cysts to pass but retains larger diameter particulate matter on the top of the screen. Most of the macroscopic fecal debris will be discarded with the funnel/filtration unit. With dense fecal samples, the flow rate will be much slower. However, do not force the sample through the filtration device by any method. Such action could damage the filtration unit requiring the sample to have to be retested. Also, scraping or forcing the sample material through the funnel/filtration device may force material through the device which will result in non-standardized sized particulates in the final sediment. This larger size material will make examination more difficult and could make coverslipping difficult.
3. The Trichrome staining method is not recommended for the examination of *Cryptosporidium spp.* The Cryptosporidium Stain Set (#0003357) is recommended.

BIBLIOGRAPHY

1. Burrows, R.B. 1965. "Microscopic Diagnosis of the Parasites of Man." Yale University Press, New Haven.
2. Garcia, L.S., R. Shimizu. 1981. "Comparison of Clinical Results for the use of Ethyl Acetate and Diethyl Ether in the Formalin-ether Sedimentation Technique Performed on Polyvinyl Alcohol Preserved Specimens." J. Clin. Microbiol. 13:709-713.
3. Garcia, L.S., M. Voge. 1990. Diagnostic Clinical parasitology. Am J. Med. Technol. 46:459-467.
4. Melvin, D.M., M.M. Brooke, 1980. Laboratory Proc. for the Diagnosis of Intestinal Parasites. U.S.D.H.E.W. 80:8282 CDC Atlanta, GA.
5. Price, D.L. 1994. Procedure Manual for the Diagnosis of Intestinal Parasites. CRC Press.
6. Ritchie, L.S. 1948. An Ether Sedimentation Technique for Routine Stool Examination. Bull. U.S. Army Med. Dept. 8:326.
7. Yang, J. and T.H. Scholten. 1977. "A Fixative for Intestinal Parasites Permitting the use of Concentration and Permanent Staining Procedures." Am. J. Clin. Pathol. 67:300-304.
8. Young, K.H., S.L. Bullock, D.M. Melvin, and C.L. Spruill. 1979 "Ethyl Acetate as a Substitute for Diethyl Ether in the Formalin-ether Sedimentation Technique." J. Clin. Microbiol. 10:852-853.

CONTACT

CalibreScientific US, Inc. offers a complete line of reagents, stains, and QC1™ Quality Control Slides for AFB, Parasitology, Bacteriology, and Mycology processing, as well as O&P collection systems and concentration devices for Parasitology. For Technical Assistance, email Technical@AlphaTecSystems.com, and for Customer Service, email Sales@AlphaTecSystems.com, or call either [+1] 800.221.6058 (USA) or [+1] 360.260.2779 between 8AM and 4PM Monday through Friday, Pacific Time.

WARRANTY

This product is warranted by CalibreScientific US, Inc. to perform as described in the labeling and literature supplied. CalibreScientific US, Inc. disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall CalibreScientific US, Inc. be liable for any consequential damages arising out of aforesaid express warranty.

TRADEMARKS:

CONSED®, CELL-BOND®, PROTO-FIX®, PRS™, and QC1™ are trademarks of CalibreScientific AMER IP LLC., 6201 Trust Dr, Holland, OH 43528.

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BIBLIOGRAFIA

1. Burrows, R.B. 1965. "Microscopic Diagnosis of the Parasites of Man." Yale University Press, New Haven.
2. Garcia, L.S., R. Shimizu. 1981. "Comparison of Clinical Results for the use of Ethyl Acetate and Diethyl Ether in the Formalin-ether Sedimentation Technique Performed on Polyvinyl Alcohol Preserved Specimens." J. Clin. Microbiol. 13:709-713.
3. Garcia, L.S., M. Voge. 1990. Diagnostic Clinical parasitology. Am J. Med. Technol. 46:459-467.
4. Melvin, D.M., M.M. Brooke, 1980. Laboratory Proc. for the Diagnosis of Intestinal Parasites. U.S.D.H.E.W. 80:8282 CDC Atlanta, GA.
5. Price, D.L. 1994. Procedure Manual for the Diagnosis of Intestinal Parasites. CRC Press.
6. Ritchie, L.S. 1948. An Ether Sedimentation Technique for Routine Stool Examinations. Bull. U.S. Army Med. Dept. 8:326.
7. Yang, J. and T.H. Scholten. 1977. "A Fixative for Intestinal Parasites Permitting the use of Concentration and Permanent Staining Procedures." Am. J. Clin. Pathol. 67:300-304.
8. Young, K.H., S.L. Bullock, D.M. Melvin, and C.L. Spruill. 1979. "Ethyl Acetate as a Substitute for Diethyl Ether in the Formalin-ether Sedimentation Technique." J. Clin. Microbiol. 10:852-853.

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MARCHI:

CONSED®, CELL-BOND®, PROTO-FIX®, PRS™, e QC1™ Slides sonomarchi di CalibreScientific AMER IP LLC., 6201 Trust Dr, Holland, OH 43528.

PRODUCT CODES

0003344 Ethyl Acetate, 1 x 500 ml
0003358 Triton X-100, 20%, 2 x 30 ml
0004043 Parasite Recovery System, 100 Sets
0004044 Parasite Recovery System, 100 Funnels
0004045 Parasite Recovery System, 100 Sets

Manufactured by CalibreScientific US, Inc.
1311 SE Cardinal Court, Suite 170
Vancouver, WA 98683 USA

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MDSS GmbH