

**THIS INFORMATION IS TO BE USED IN CONJUNCTION  
WITH DFU L004302, NAC-PAC RED****PROCESSING BLOODY AFB SPECIMENS**

If an AFB specimen submitted for processing contains a significant amount of blood it may be difficult to visualize the color change from pink to colorless when adding the AFB Neutralization Buffer (XPR-PLUS® or NPC-67®) to the specimen and NAC-PAC RED solution after digestion. This is a modification to ensure proper neutralization when color change cannot be visualized. When adding the NAC-PAC RED to the bloody specimen, the NaOH in the solution will lyse the red blood cells and the residual hemoglobin can mask the visualization of the pH indicator in the solution in the neutralization step.

1. To the sterile 50 ml centrifuge tube, add the NAC-PAC RED/NALC solution in the following amounts: for specimens 1-5 ml, add a volume equal to that of the specimen volume; to specimens 6-7 ml, add 5 ml of base solution; for specimens 8-10 ml, add an equal volume of base digestant and split the specimen after step 3 equally into two centrifuge tubes, proceed with steps 4-6 and then combine the sediments from both tubes into one centrifuge tube and proceed with step 8.
2. Mix each specimen on a vortex until liquefied.
3. Allow each specimen to stand for 15 minutes, vortexing every 5 minutes during this step. Each specimen should stand for 15 minutes, but no more than 20 minutes.
4. Add the Neutralization Buffer up to the 50mL mark on the centrifuge tube to bring the solution pH below 8.1. This will ensure optimal recovery of any Mycobacterium that may be present in the specimen.
5. Centrifuge the specimen tubes at 3000 xg for 15 minutes. It is recommended but not required to use a refrigerated centrifuge.
6. Working in a biosafety hood, pour off all of the supernatant into a splash-proof container holding an appropriate disinfectant.
7. Resuspend the pellet with 0.5-1.0 ml of PRB™ Pellet Resuspension Buffer.
8. Prepare appropriate smears and/or inoculate culture media according to laboratory protocols.

The lysed red blood cells will also result in a significant amount of cell debris in the solution. The debris may attach to the sides of the centrifuge tube due to the natural charge associated with the plastic.

**PROCESSING SMALL VOLUME AFB SPECIMENS**

Small volume specimens with corresponding low post-neutralization volumes can make centrifuge balancing difficult. This is a modification to accommodate small volumes (under 5 mL).

1. Add **sterile** saline to the sample in a sterile 50 ml centrifuge tube to reach a total volume of 5 mL.
2. Add 5 mL of NAC-PAC RED containing N-acetyl-L cysteine.
3. Mix each specimen on a vortex until liquefied.
4. Allow each specimen to stand for 15 minutes, vortexing every 5 minutes during this step. Each specimen should stand for 15 minutes, but no more than 20 minutes.
5. Centrifuge the specimen tubes at 3000 xg for 15 minutes. It is recommended but not required to use a refrigerated centrifuge.
6. Working in a biosafety hood, pour off all the supernatant into a splash-proof container holding an appropriate disinfectant.
7. Resuspend the pellet with 0.5-1.0 ml of PRB.
8. Prepare appropriate smears and/or inoculate culture media according to laboratory protocols.

**PROCESSING SPECIMENS CONTAMINATED WITH  
PSUEDOMONAS SPP.**

Specimens contaminated with *Pseudomonas* spp. will need additional treatment with 5% Oxalic Acid to prevent overgrowth of *Pseudomonas* spp. Refer to the 5% Oxalic Acid Directions For Use (DFU L003447) for complete instructions. For additional information on the pH effects of the oxalic acid procedure and the appropriate buffering requirements, call Alpha-Tec Technical Services or the Alpha-Tec Sales Department.

**RE-PROCESSING CONTAMINATED AFB SPECIMENS:**

If there is non-AFB overgrowth on the test media or rapid growth detection systems, the specimens must be re-processed to remove the contaminating bacteria.

1. Add **sterile** saline or Pellet Resuspension Buffer to the sample in a sterile 50 ml centrifuge tube to reach a total volume of 5 mL.
2. Add 5 mL of NAC-PAC RED containing N-acetyl-L cysteine.
3. Mix each specimen on a vortex until liquefied.
4. Allow each specimen to stand for 15 minutes, vortexing every 5 minutes during this step. Each specimen should stand for 15 minutes and can be increased up to 20 minutes to ensure adequate decontamination.
5. Centrifuge the specimen tubes at 3000 xg for 15 minutes. It is recommended but not required to use a refrigerated centrifuge.
6. Working in a biosafety hood, pour off all of the supernatant into a splash-proof container holding an appropriate disinfectant.
7. Resuspend the pellet with 0.5-1.0 ml of PRB.
8. Prepare appropriate smears and/or inoculate culture media according to laboratory protocols.

**PROCESSING STOOL AFB SPECIMENS:**

If a stool specimen is submitted to the laboratory for AFB processing it may be difficult to visualize the color change from pink to colorless when adding the AFB Neutralization Buffer (XPR-PLUS or NPC-67) due to the opacity of fecal material. This is a modification to the digestion and decontamination procedure (NAC-PAC RED DFU L004302) to accommodate stool specimens.

1. In a sterile 50 ml centrifuge tube, suspend approximately 1 gram of fecal material in 5 ml of saline. Vortex to mix well.
2. Add 5 ml of NAC-PAC RED containing N-acetyl-L-cysteine and allow the specimen to liquefy and decontaminate for 15 minutes, vortexing every 5 minutes.
3. Add enough Neutralization Buffer to neutralize the NAC-PAC RED (liquid goes from red to colorless). If you are unable to see the color change due to the fecal material, add the Neutralization Buffer up to the 50 ml mark on the centrifuge tube to bring the pH below 8.1.
4. Centrifuge the specimen tube at 3000 xg for 15 minutes. It is recommended, but not required to use a refrigerated centrifuge.
5. Working in a biosafety hood, pour off all of the supernatant into a splash-proof container holding an appropriate disinfectant.
6. Resuspend the pellet with 0.5 – 1.0 ml of Pellet Resuspension Buffer.
7. Prepare appropriate smears and/or inoculate culture media according to laboratory protocols.
8. Save the remaining resuspended pellet in case of excessive breakthrough contamination. If the breakthrough contamination is unacceptable, re-process the pellet following the directions in the section “Re-processing Contaminated Specimens” in this document, increasing the decontamination time to 20 minutes.
9. Due to the inherent high bioburden of stool specimens, you may need to use a higher (up to 4.0%) concentration of NAC-PAC RED solution to help control breakthrough contamination.



[+1] 800.221.6058 United States  
[+1] 360.260.2779 International

Directions For Use for the following:  
(Supplement Information)

**NAC-PAC<sup>®</sup> RED**

---

**CONTACT**

Alpha-Tec Systems, Inc. offers a complete line of reagents, stains, and QC1™ Quality Control Slides for AFB, Parasitology, Bacteriology, and Mycology processing, as well as O&P collection systems and concentration devices for Parasitology. For Technical Assistance, email [Technical@AlphaTecSystems.com](mailto:Technical@AlphaTecSystems.com), and for Customer Service, email [Sales@AlphaTecSystems.com](mailto:Sales@AlphaTecSystems.com), or call either [+1] 800.221.6058 (USA) or [+1] 360.260.2779 between 8AM and 4PM Monday through Friday, Pacific Time.

**WARRANTY**

This product is warranted by Alpha-Tec Systems, Inc. to perform as described in the labeling and literature supplied. Alpha-Tec Systems, Inc. disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall Alpha-Tec Systems, Inc. be liable for any consequential damages arising out of aforesaid express warranty.

**TRADEMARKS**

NAC-PAC RED<sup>®</sup>, NPC-67<sup>®</sup>, PRB™ and XPR-PLUS<sup>®</sup> are trademarks of Alpha-Tec Systems, Inc., 1311 SE Cardinal Court, Suite 170, Vancouver, WA 98683 USA.