**TABLE OF CONTENTS**

English...	1	German...	7
French...	3	Product Codes...	9
Italian...	5	Glossary of Symbols...	9

**INTENDED USE**

NAC-PAC® UNITARY™ Individualized Specimen Reagents for Mycobacteria Processing contains the necessary reagents [NAC-PAC RED, NPC-67® Neutralizing Buffer, and PRB™ Pellet Resuspension Buffer] for use in the qualitative digestion and decontamination procedure of clinical specimens for the recovery of *Mycobacterium* spp. UNITARY Specimen Tubes are 50 ml centrifuge tubes containing a N-acetyl-L-cysteine (NALC) tablet and are used in conjunction with the NAC-PAC UNITARY system.

**SUMMARY**

The decontamination and digestion procedure, utilizing the compound NALC combined with sodium hydroxide and sodium citrate solution, results in increased yields of tubercle bacilli. The NALC procedure utilizes N-acetyl-L-cysteine as a mucolytic compound by disrupting chemical bonds in mucus. The sodium hydroxide acts as a bacterial decontaminant and the sodium citrate stabilizes the NALC by chelating (binding) any heavy metal ions. Since the sodium hydroxide has a pH of approximately 13.00, it will kill bacteria (including mycobacteria) after 15–20 minutes of exposure. Therefore, timing of the decontamination is critical to limit the amount of *Mycobacterium* spp. killed by the basic pH. A pH indicator is incorporated into the NAC-PAC RED decontamination reagent to monitor the pH throughout the decontamination and buffering procedure. Bringing the pH to a neutral range stops the decontamination process. The NPC-67 is used to neutralize the NAC-PAC RED following the appropriate digestion and decontamination time, resulting in a pH of below 8.10. Studies have documented that pH values above 8.10 are toxic to *Mycobacterium* spp., including *Mycobacterium tuberculosis*. Following the decanting step, PRB is added to achieve a tight neutral pH value (6.80–7.10) in the specimen sediment, optimizing mycobacteria recovery.

**FOR IN VITRO DIAGNOSTIC USE ONLY**
**PRECAUTIONS**

NAC-PAC RED contains a caustic chemical (sodium hydroxide). Use appropriate care in the handling of this reagent. All clinical specimens submitted for the diagnosis of tuberculosis and other *Mycobacterium* spp. must be treated with appropriate care so as not to contaminate other specimens or laboratory personnel. Use all approved and regulated equipment for processing and detection procedures.

**STABILITY AND STORAGE**

Reagents in NAC-PAC UNITARY are stable to the stated expiration date when stored at the required temperature. Store at room temperature (15–30°C). Do not freeze or heat above 30°C.

**USER QUALITY CONTROL**

Any product showing cloudiness, turbidity, precipitation, or discoloration should be discarded. Quality controlled microorganisms should be utilized to verify procedures, media and reagents as appropriate for your laboratory's applicable regulatory agency or local procedural guidelines.

**SPECIMEN COLLECTION AND PREPARATION**

Appropriate specimens for the detection of *Mycobacterium* spp. should be collected according to prescribed standards and delivered to the laboratory in a safe and timely manner. Refer to local procedural guidelines for this information.

**PROCEDURE**

**Materials Provided with NAC-PAC UNITARY (#0004810):** Contains reagents for the processing of 20 patient specimens. There are twenty

bottles of each reagent: NAC-PAC RED, NPC-67 Neutralizing Buffer, and PRB Pellet Resuspension Buffer.

**Materials Provided with UNITARY Specimen Tubes (#0003928):** Twenty 50 ml centrifuge tubes, each containing a NALC tablet.

**Materials Not Provided:** Centrifuge, vortex mixer, sterile pipettes, microscope slides, TB media, Alpha-Tec CELL-BOND® Slides.

**SPECIMEN PROCESSING**

1. Work in sets equivalent to your centrifuge capacity. In a biosafety hood, remove cap from the UNITARY Specimen Tube containing a NALC tablet and place up to 7 ml of patient sample inside. If the sample is larger than 7 ml, split the sample into multiple tubes. If the sample is refrigerated, allow the sample to come to room temperature in the UNITARY Specimen Tube.
2. Once the sample is at room temperature, add the contents of NAC-PAC RED as follows:
  - a. For specimens 1-4 ml, add an equal volume of NAC-PAC RED.
  - b. For specimens 5-7 ml, add the entire contents of NAC-PAC RED.
3. Replace and tighten cap on UNITARY Tube, and vortex for 30 seconds.
4. Allow each specimen to stand for 15 minutes, but no longer than 20 minutes. Vortex every 5 minutes during this step.
5. Add NPC-67 to the 50 ml mark of the tube. Effective neutralization is indicated by a color change from red/pink to colorless. Replace and tighten cap, and swirl by hand or invert gently 3–5 times to mix.
6. Centrifuge at 3000 xg for 15 minutes. It is recommended, but not required, to use a refrigerated centrifuge. Each laboratory must check the centrifuge head radius, and use an appropriate nomogram for proper speed selection (rpm) to achieve the desired relative centrifugal field of 3000 xg.
7. Working in a biosafety hood, pour off all supernatant into a splash-proof container holding an appropriate disinfectant. Use an appropriate disinfectant to disinfect any contamination on the lip of the UNITARY Tube. Do not allow the disinfectant to run down the inside of the tube.
8. Resuspend the pellet with 1 ml of PRB.
9. Mix the sediment and buffer well and make smears for acid-fast staining. Use adhesive CELL-BOND Slides or appropriate sterile albumin adhesive solutions to attach the specimen to the slide. Dry the smears and proceed with acid-fast staining per the manufacturer's directions. **NOTE:** An acid-fast stain control slide should be stained in conjunction with the patient smears to verify the staining technique and components (#0003240, #0003245, or #0003247 QC1™ AFB Slide).
10. Place two drops of the sediment onto the surface of each of the TB media used and inoculate to rapid broth detection systems, per the manufacturer's recommendations. **NOTE:** A BAP or TSA contamination control plate can be inoculated and incubated at 35–37°C for 48 hrs.
11. To the unused portion of the specimen, add the remainder of PRB. Save the balance of the pellet specimen and refrigerate at 2–8°C for further diagnostic procedures or reprocessing, if necessary.

**PROCEDURE NOTES**

1. The components of NAC-PAC UNITARY and UNITARY Specimen Tubes have been validated for use with multiple molecular diagnostic methods and systems. For more information regarding compatibility with specific methods or systems, contact Alpha-Tec Technical Services.
2. Specimens contaminated with *Pseudomonas* spp. will need additional treatment with Oxalic Acid. Refer to the Oxalic Acid Directions For Use for complete instructions, or contact Alpha-Tec Technical Services or the Alpha-Tec Sales Department for information on the pH effects of the oxalic acid procedure and the appropriate buffering requirements.
3. Bloody Specimens: Following the decontamination of the specimen with NAC-PAC RED, bloody specimens may remain pink after the addition of the NPC-67 Neutralizing Buffer due to the residual hemoglobin in the specimen. Adding NPC-67 to the 50 ml mark will

ensure complete neutralization. For additional information, contact Alpha-Tec Technical Services.

#### EXPECTED RESULTS

If *Mycobacterium* spp. are present in the clinical specimen and processed according to the procedures listed within this document, the recovery of cultivable, viable and clinically significant *Mycobacterium* spp. can be expected.

#### LIMITATION OF PROCEDURES

Timing of the decontamination step, proper buffering, speed and timing of the centrifugation step, proper decanting and the addition of the PRB to the pellet are vital to the recovery of *Mycobacterium* spp. Failure to follow the listed procedures may result in decreased numbers or total loss of *Mycobacterium* spp. resulting in an inaccurate culture report.

#### SPECIFIC PERFORMANCE CHARACTERISTICS

The products within NAC-PAC UNITARY were tested on clinical samples and recovered all culture-appropriate *Mycobacterium* spp. when the designated procedures were followed.

#### BIBLIOGRAPHY

1. Kent, P., Kubica, G.P., *Public Health Mycobacteriology: A Guide for the Level III Laboratory*. Centers for Disease Control and Prevention (CDC). 1985.
2. Kubica, G.P., et al., "Sputum Digestion and Decontamination with N-acetyl-L-cysteine-Sodium Hydroxide for Culture of Mycobacteria." *Am Rev Respir Dis*. 1963. May 87:775-779.
3. Lennette, E.H., et al., *Manual of Clinical Microbiology Third Addition*. American Society for Microbiology. 1980.
4. Murray, P., *Manual of Clinical Microbiology Eighth Addition*. American Society for Microbiology. 2003.
5. *Mycobacterium tuberculosis: Assessing Your Laboratory*. Centers for Disease Control and Prevention (CDC). March, 1995.
6. Vestal, A.L., *Procedures for the Isolation and Identification of Mycobacteria*. Centers for Disease Control and Prevention (CDC), Atlanta, GA. 1975.
7. Yegian, D., Budd, V. "Toxic Effect of Sodium Hydroxide on Tubercle Bacilli." *Am. J. Clinical Pathology*. 1952. 22:456:460.

#### CONTACT

For technical assistance email [Technical@AlphaTecSystems.com](mailto:Technical@AlphaTecSystems.com) and for Customer Service, email [Sales@AlphaTecSystems.com](mailto:Sales@AlphaTecSystems.com) or call [+1] 800.221.6058 or [+1] 360.260.2779 between 8 am and 4 pm Monday through Friday, Pacific Time.

#### WARRANTY

Alpha-Tec Systems, Inc. warrants this product to perform as described in the labeling and literature supplied. Alpha-Tec Systems, Inc. disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall Alpha-Tec Systems, Inc. be liable for any consequential damages arising out of aforesaid express warranty.

#### TRADEMARKS

CELL-BOND®, NAC-PAC®, NPC-67®, PRB™, QC1™, and UNITARY™ are trademarks of Alpha-Tec Systems, Inc., 1311 SE Cardinal Court, Suite 170, Vancouver, WA 98683 USA.

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## BIBLIOGRAPHIE

1. Kent, P., Kubica, G.P., *Public Health Mycobacteriology: A Guide for the Level III Laboratory*. Centers for Disease Control and Prevention (CDC). 1985.
2. Kubica, G.P., et al., "Sputum Digestion and Decontamination with N-acetyl-L-cysteine-Sodium Hydroxide for culture of Mycobacteria." *Am Rev Respir Dis*. 1963. May 87:775-779.
3. Lennette, E.H., et al., *Manual of Clinical Microbiology Third Addition*. American Society for Microbiology. 1980.
4. Murray, P., *Manual of Clinical Microbiology Eighth Addition*. American Society for Microbiology. 2003.
5. *Mycobacterium tuberculosis: Assessing Your Laboratory*. Centers for Disease Control and Prevention (CDC). March, 1995.
6. Vestal, A.L., *Procedures for the Isolation and Identification of Mycobacteria*. Centers for Disease Control and Prevention (CDC), Atlanta, GA. 1975.
7. Yegian, D., Budd, V. "Toxic Effect of Sodium Hydroxide on Tubercle Bacilli." *Am. J. Clinical Pathology*. 1952. 22:456:460.

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**BIBLIOGRAFIA**

1. Kent, P., Kubica, G.P., *Public Health Mycobacteriology: A Guide for the Level III Laboratory*. Centers for Disease Control and Prevention (CDC). 1985.
2. Kubica, G.P., et al., "Sputum Digestion and Decontamination with N-acetyl-L-cysteine-Sodium Hydroxide for Culture of Mycobacteria." *Am Rev Respir Dis*. 1963. May 87:775-779.
3. Lennette, E.H., et al., *Manual of Clinical Microbiology Third Edition*. American Society for Microbiology. 1980.
4. Murray, P., *Manual of Clinical Microbiology Eighth Edition*. American Society for Microbiology. 2003.
5. *Mycobacterium tuberculosis: Assessing Your Laboratory*. Centers for Disease Control and Prevention (CDC). March, 1995.
6. Vestal, A.L., *Procedures for the Isolation and Identification of Mycobacteria*. Centers for Disease Control and Prevention (CDC), Atlanta, GA. 1975.
7. Yegian, D., Budd, V. "Toxic Effect of Sodium Hydroxide on Tubercle Bacilli." *Am. J. Clinical Pathology*. 1952. 22:456:460.

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#### LITERATUR

- Kent, P., Kubica, G.P., *Public Health Mycobacteriology: A Guide for the Level III Laboratory*. Centers for Disease Control and Prevention (CDC). 1985.
- Kubica, G.P., et al., "Sputum Digestion and Decontamination with N-acetyl-L-cysteine-Sodium Hydroxide for Culture of Mycobacteria." *Am Rev Respir Dis*. 1963. May 87:775-779.
- Lennette, E.H., et al., *Manual of Clinical Microbiology Third Edition*. American Society for Microbiology. 1980.
- Murray, P., *Manual of Clinical Microbiology Eighth Edition*. American Society for Microbiology. 2003.
- Mycobacterium tuberculosis: Assessing Your Laboratory*. Centers for Disease Control and Prevention (CDC). March, 1995.

- Vestal, A.L., *Procedures for the Isolation and Identification of Mycobacteria*. Centers for Disease Control and Prevention (CDC), Atlanta, GA. 1975.
- Yegjian, D., Budd, V. "Toxic Effect of Sodium Hydroxide on Tubercle Bacilli." *Am. J. Clinical Pathology*. 1952. 22:456:460.

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**PRODUCT CODES**

0003928 NAC-PAC UNITARY Specimen Tubes, 20 tubes  
0004810 NAC-PAC UNITARY Reagents (3.0%), 20 tests

Manufactured by Alpha-Tec Systems, Inc.  
1311 SE Cardinal Court, Suite 170  
Vancouver, WA 98683 USA

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MDSS GmbH