TABLE OF CONTENTS

English...	1	German...	9
French...	3	Product Codes...	11
Spanish...	5	Glossary of Symbols...	11
Italian...	7		

INTENDED USE

NAC-PAC® EA3 contains the necessary reagents—NAC-PAC RED, NALC tablets or powder, NPC-67® Neutralizing Buffer, and PRB™ Pellet Resuspension Buffer—for use in the qualitative digestion and decontamination procedure of clinical specimens for the recovery of *Mycobacterium* spp.

SUMMARY

The decontamination and digestion procedure, utilizing the compound N-acetyl-L-cysteine (NALC) combined with sodium hydroxide and sodium citrate solution, results in increased yields of tubercle bacilli. The NALC procedure utilizes N-acetyl-L-cysteine as a mucolytic compound by disrupting chemical bonds in mucus. The sodium hydroxide acts as a bacterial decontaminant and the sodium citrate stabilizes the NALC by chelating (binding) any heavy metal ions. Since the sodium hydroxide has a pH of approximately 13.00, it will kill bacteria (including mycobacteria after 15–20 minutes of exposure). Therefore, timing of the decontamination is critical to limit the amount of *Mycobacterium* spp. killed by the basic pH. A pH indicator is incorporated into the NAC-PAC RED decontamination reagent to monitor the pH throughout the decontamination and buffering procedure. Bringing the pH to a neutral range stops the decontamination process. The NPC-67 is used to neutralize the NAC-PAC RED following the appropriate digestion and decontamination time, resulting in a pH below 8.10. Studies have documented that pH values above 8.10 are toxic to *Mycobacterium* spp., including *Mycobacterium tuberculosis*. Following the decanting step, PRB is added to achieve a tight neutral pH value (6.80–7.10) in the specimen sediment, optimizing mycobacteria recovery.

FOR IN VITRO DIAGNOSTIC USE ONLY

PRECAUTIONS

NAC-PAC RED contains a caustic chemical (sodium hydroxide). Use appropriate care in the handling of this reagent. All clinical specimens submitted for the diagnosis of tuberculosis and other *Mycobacterium* spp. must be treated with appropriate care to avoid contaminating other specimens or laboratory personnel. Use approved and regulated equipment for processing and detection procedures.

STABILITY AND STORAGE

Reagents in the NAC-PAC EA3 are stable to the stated expiration date when stored at the required temperature. Store at room temperature (15–30°C).

QUALITY CONTROL

Any reagent showing cloudiness, turbidity, precipitation, or discoloration should be discarded. Quality controlled microorganisms should be utilized to verify procedures, media, and reagents as appropriate for your laboratory's applicable regulatory agency or local procedural guidelines.

SPECIMEN COLLECTION AND PREPARATION

Appropriate specimens for the detection of *Mycobacterium* spp. should be collected according to prescribed standards and delivered to the laboratory in a safe and timely manner. Refer to local procedural guidelines for this information.

PROCEDURE

Materials Provided: NAC-PAC RED, NALC tablets or powder, NPC-67 Neutralizing Buffer, PRB Pellet Resuspension Buffer.

Materials Not Provided: Centrifuge and 50 ml centrifuge tubes, vortex mixer, sterile pipettes, microscope slides, TB growth media, CELL-BOND® Slides.

SPECIMEN PROCESSING

1. Prior to opening the specimen tube(s), line up specimens (collected in centrifuge tubes) in a rack inside an appropriate biosafety hood.
2. Working in sets equivalent to your centrifuge capacity, loosen specimen container caps. To prevent the possibility of cross-contamination, work with one specimen at a time.
3. Open the bottle labeled "NALC" and place one tablet in each centrifuge tube containing a specimen. Or, add the NALC powder to the bottle labeled "NAC-PAC RED," tightly cap the bottle, and shake well or vortex to dissolve the NALC. **Once dissolved, the NAC-PAC RED/NALC solution will be stable for 72 hours. Store any unused portion at 2–8°C for up to 72 hours. Allow the refrigerated portion to come to room temperature prior to use.**
4. To the specimen tube containing the NALC tablet, add NAC-PAC RED (or the NAC-PAC RED/NALC solution) as follows:
 - a. For specimens **1–5 ml add a volume of NAC-PAC RED equal to that of the specimen volume.**
 - b. For specimens **6–7 ml add 5 ml of NAC-PAC RED.**
 - c. For specimens **8–10 ml add an equal volume of NAC-PAC RED and then split the specimen after step 6** equally into two centrifuge tubes, proceed with steps 7–9 and then combine the sediments from both tubes into one centrifuge tube and proceed with step #10. Following this protocol will help achieve effective decontamination while also allowing for proper neutralization.
5. Tighten the caps on the centrifuge tubes. Mix each specimen on a vortex until liquefied and the NALC tablet has dissolved (30 seconds per specimen).
6. Allow each specimen to stand for 15–20 minutes. Vortex every five minutes during this step.
7. Fill each tube with NPC-67 until effective neutralization is indicated by a color change from red/pink to colorless. Once a colorless point has been reached, do not continue to add NPC-67 to the sample. Tighten cap and swirl by hand to mix. **NOTE:** Use each bottle of NPC-67 on a single specimen only. Using one bottle for multiple specimens can lead to cross-contamination and potentially erroneous results.
8. Centrifuge the specimen tubes at 3000 xg for 15 minutes. It is recommended but not required to use a refrigerated centrifuge. Each laboratory must check the centrifuge head radius, and use an appropriate nomogram for proper speed selection [rpm] to achieve the desired relative centrifugal field of 3000 xg.
9. Working in a biosafety hood, pour off all supernatant into a splash-proof container holding an appropriate disinfectant. Use an appropriate disinfectant to disinfect any contamination on the lip of the specimen tube. Do not allow the disinfectant to run down inside the specimen tube.
10. Resuspend the pellet with 0.5–1.0 ml of PRB. Do not resuspend the pellet with NPC-67, water, or saline. **NOTE:** To maximize time to detection for rapid growth automated detection systems, resuspend the pellet with 1.0 ml of PRB. Depending on the needs of your laboratory, the pellet may be resuspended with 0.5 ml of PRB to create a more concentrated sample for increased acid-fast smear staining sensitivity. Once the smears have been made, add an additional 1.0 ml of PRB to inoculate rapid broth detection systems and other media.
11. Mix the sediment and buffer well and inoculate the liquid broth for your automated detection equipment per the manufacturer's instructions.
12. Place two drops of the sediment onto the surface of each of the TB media used. **NOTE:** A contamination control plate [BAP or TSA] can be inoculated at this point and incubated at 35–37°C for 48 hours.
13. Make smears for acid-fast staining. Use adhesive CELL-BOND slides or sterile albumin adhesive solutions to attach the specimen to the slide. Dry the smears and proceed with appropriate acid-fast staining per the manufacturer's directions. **NOTE:** An acid-fast stain control slide should be stained in conjunction with the patient smears to verify the staining technique and components. Contact Alpha-Tec for a complete list of acid-fast stains and control slides.

14. To the unused portion of the specimen, add the balance of the PRB and refrigerate at 2–8°C to save for further diagnostic procedures or reprocessing, if necessary.

PROCEDURE NOTES

1. Alpha-Tec NAC-PAC EA3 has been validated for use with multiple molecular diagnostic methods and systems. For more information regarding compatibility with specific methods or systems, contact Alpha-Tec Technical Services.
2. Small volume specimens: Small volume specimens with correspondingly low post-neutralization volumes can make centrifuge balancing difficult. If your laboratory frequently encounters small volume specimens, it is acceptable to add sterile saline to the sample to reach a combined volume of 5 ml prior to the addition of the NALC tablet. In this case, the sample should be decontaminated with 5 ml of NAC-PAC RED solution. This will increase the final post-neutralization specimen volume making centrifuge balancing easier.
3. Specimens contaminated with *Pseudomonas* spp. will need additional treatment with 5% Oxalic Acid (OxA® Kit #0004805). Refer to the OxA Kit Directions For Use for complete instructions, or call Alpha-Tec Technical Services or the Alpha-Tec Sales Department for information on the pH effects of the oxalic acid procedure and the appropriate buffering requirements.
4. Bloody Specimens: Following the decontamination of the specimen with NAC-PAC RED, bloody specimens may remain pink after the addition of the NPC-67 due to the residual hemoglobin in the specimen. If the color change cannot be visualized due to hemoglobin, add the NPC-67 up to the 50 ml mark to ensure complete neutralization. For additional information, contact Alpha-Tec Technical Services.

EXPECTED RESULTS

If *Mycobacterium* spp. are present in the clinical specimen and processed according to the procedures listed within this document, the recovery of cultivable, viable, and clinically significant *Mycobacterium* spp. can be expected.

LIMITATIONS OF PROCEDURES

Timing of the decontamination step, proper buffering, speed and timing of the centrifugation step, proper decanting, and the addition of the PRB to the pellet are vital to the recovery of *Mycobacterium* spp. Failure to follow the listed procedures may result in decreased numbers or total loss of *Mycobacterium* spp. resulting in an inaccurate culture report.

SPECIFIC PERFORMANCE CHARACTERISTICS

The products within NAC-PAC EA3 were tested on clinical samples and recovered all culture-appropriate *Mycobacterium* spp. when the designated procedures were followed.

BIBLIOGRAPHY

1. Kent, P, Kubica, GP. Public Health Mycobacteriology: A Guide for the Level III Laboratory. Centers for Disease Control and Prevention (CDC). 1985.
2. Kubica, GP, et al. Sputum digestion and decontamination with N-acetyl-L-cysteine-Sodium Hydroxide for culture of mycobacteria. Am Rev Respir Dis. 1963. 87:775–779.
3. Kubica, GP, et al. Comments on the use of the new mucolytic agent, N-acetyl-L-cysteine, as a sputum digestant for the isolation of mycobacteria. Am Rev Respir Dis. 1964. 89:284–286.
4. Lennette, EH, et al. Manual of Clinical Microbiology. American Society for Microbiology. 3rd ed. 1980.
5. Murray, P, Manual of Clinical Microbiology. American Society for Microbiology. 1985.
6. Mycobacterium Tuberculosis: Assessing Your Laboratory. Centers for Disease Control and Prevention (CDC). March, 1995.
7. Vestal, AL. Procedures for the Isolation and Identification of Mycobacteria. CDC. Atlanta, GA. 1975.
8. Yegian, D, Budd V. Toxic effect of sodium hydroxide on tubercle bacilli. Am J Clin Pathol 1952. 22:456–460.

CONTACT

For Technical Assistance email Technical@AlphaTecSystems.com and for Customer Service, email Sales@AlphaTecSystems.com or call [+1] 800.221.6058 or [+1] 360.260.2779 between 8 am and 4 pm Monday through Friday, Pacific Time.

WARRANTY

Alpha-Tec Systems, Inc. warrants this product to perform as described in the labeling and literature supplied. Alpha-Tec Systems, Inc. disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall Alpha-Tec Systems, Inc. be liable for any consequential damages arising out of aforesaid express warranty.

TRADEMARKS

CELL-BOND®, NAC-PAC®, NPC-67®, OxA®, and PRB™ are trademarks of Alpha-Tec Systems, Inc., 1311 SE Cardinal Court, Suite 170, Vancouver, WA 98683 USA.

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BIBLIOGRAFIA

- Kent, P, Kubica, GP. Public Health Mycobacteriology: A Guide for the Level III Laboratory. Centers for Disease Control and Prevention (CDC). 1985.
- Kubica, GP, et al. Sputum digestion and decontamination with N-acetyl-L-cysteine-Sodium Hydroxide for culture of mycobacteria. Am Rev Respir Dis. 1963. 87:775–779.
- Kubica, GP, et al. Comments on the use of the new mucolytic agent, N-acetyl-L-cysteine, as a sputum digestant for the isolation of mycobacteria. Am Rev Respir Dis. 1964. 89:284–286.
- Lennette, EH, et al. Manual of Clinical Microbiology. American Society for Microbiology. 3rd ed. 1980.
- Murray, P, Manual of Clinical Microbiology. American Society for Microbiology. 1985.
- Mycobacterium Tuberculosis: Assessing Your Laboratory. Centers for Disease Control and Prevention (CDC). March, 1995.
- Vestal, AL. Procedures for the Isolation and Identification of Mycobacteria. CDC. Atlanta, GA. 1975.
- Yegian, D, Budd V. Toxic effect of sodium hydroxide on tubercle bacilli. Am J Clin Pathol 1952. 22:456–460.

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Manufactured by Alpha-Tec Systems, Inc.
MDSS GmbH